


# TopSpin

- Advanced NMR Experiments  
User Manual  
Version 005



Copyright © by Bruker Corporation

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form, or by any means without the prior consent of the publisher. Product names used are trademarks or registered trademarks of their respective holders.

© November 23, 2020 Bruker Corporation

Document Number: 1000068378

P/N: H147756

---

# Contents

<b>1</b>	<b>About This Manual .....</b>	<b>9</b>
1.1	Symbols and Conventions .....	9
<b>2</b>	<b>Introduction.....</b>	<b>11</b>
2.1	General .....	11
2.2	Disclaimer .....	11
<b>3</b>	<b>The TopSpin Interface.....</b>	<b>13</b>
<b>4</b>	<b>1D Solvent Suppression Experiments.....</b>	<b>27</b>
4.1	Introduction .....	27
4.2	Samples .....	27
4.3	Preparation Experiment .....	27
4.3.1	Acquisition.....	29
4.3.2	Processing .....	29
4.4	1D Solvent Suppression with Pre-saturation .....	30
4.4.1	Parameter Setup.....	30
4.4.2	Fine Tuning .....	31
4.4.3	Acquisition.....	33
4.4.4	Processing .....	33
4.5	1D Solvent Suppression with Presaturation and Composite Pulses.....	35
4.5.1	Parameter Setup.....	35
4.5.2	Acquisition.....	35
4.5.3	Processing .....	35
4.6	1D Solvent Suppression Using the Noesy Sequence.....	36
4.6.1	Parameter Setup.....	37
4.6.2	Processing .....	37
4.7	1D Solvent Suppression with WATERGATE .....	39
4.7.1	Parameter Setup.....	39
4.7.2	Processing .....	40
4.8	1D Solvent Suppression with Excitation Sculpting.....	41
4.8.1	Parameter Setup.....	41
4.8.2	Processing .....	42
4.9	1D Solvent Suppression with WET .....	43
4.9.1	Sample.....	43
4.9.2	Preparation Experiment .....	44
4.9.3	Acquisition.....	44
4.9.4	Processing .....	44
4.9.5	Selective Excitation Region Setup .....	45
4.9.6	Processing .....	48
<b>5</b>	<b>2D Gradient Experiments.....</b>	<b>51</b>
5.1	Introduction .....	51
5.2	Setting up the COSY Experiment .....	51
5.2.1	Limit Setting .....	52

5.2.2	Acquisition.....	53
5.2.3	Processing .....	53
5.3	Setting up the COSYGP Experiment .....	55
5.3.1	Limit Setting .....	56
5.3.2	Acquisition.....	57
5.3.3	Processing .....	57
5.3.4	N- to P-Type Coherence Selection .....	57
5.4	Setting up the Multiple-Quantum Filtered COSY Experiment .....	58
5.4.1	Preparation Experiment .....	59
5.4.2	Setting up the MQF-COSY Experiment .....	60
5.4.3	Limit Setting .....	60
5.4.4	Acquisition.....	61
5.4.5	Processing .....	62
5.4.6	Plotting .....	62
5.5	Setting up the <sup>13</sup> C-HMBCGP Experiment .....	63
5.5.1	Limit Setting .....	65
5.5.2	Acquisition.....	66
5.5.3	Processing .....	66
5.6	Setting up the <sup>15</sup> N-HMBCGP from the <sup>13</sup> C-HMBCGP .....	66
5.6.1	Limit Setting .....	66
5.7	Setting up the <sup>15</sup> N-HMBC Experiment .....	68
5.7.1	Acquisition.....	70
5.7.2	Processing .....	70
5.8	2D <sup>1</sup> H / <sup>31</sup> P Gradient HMBC Experiment .....	71
5.8.1	Sample .....	71
5.8.2	<sup>1</sup> H Reference Experiment .....	71
5.8.3	<sup>31</sup> P Reference Experiment .....	72
5.8.4	Acquisition.....	73
5.8.5	Processing .....	73
5.8.6	Setting up the HMBC Experiment .....	73
5.8.7	Limit Setting .....	76
5.8.8	Acquisition.....	79
5.8.9	Processing .....	79
5.8.10	Plotting .....	80
<b>6</b>	<b>1D Experiments Using Shaped Pulses .....</b>	<b>81</b>
6.1	Introduction .....	81
6.2	Sample .....	81
6.3	1D Selective COSY Using the Flow Bar Tools .....	82
6.3.1	Introduction .....	82
6.3.2	Reference Spectrum .....	82
6.3.3	Selective Excitation Region Set Up .....	83
6.3.4	Processing .....	86
6.3.5	Plotting Two Spectra on the Same Page .....	87
6.4	1D Selective NOESY Using the Flow Bar Tools .....	88
6.4.1	Introduction .....	88
6.4.2	Reference Spectrum .....	89
6.4.3	Selective Excitation Region Set Up .....	89

6.4.4	Processing .....	92
6.4.5	Plotting Two Spectra on the Same Page .....	93
6.5	1D Selective TOCSY Using the Flow Bar Tools .....	94
6.5.1	Reference Spectrum .....	95
6.5.2	Selective Excitation Region Set Up .....	95
6.5.3	Processing .....	98
6.5.4	Plotting Two Spectra on the Same Page .....	99
6.5.5	Plotting All 4 Experiments on the Same Page .....	100
6.6	1D Selective COSY Experiment Using the On- Resonance Option .....	101
6.6.1	Reference Spectrum .....	101
6.6.2	Selective Excitation Region Set Up .....	102
6.6.3	Setting Up the Selective COSY .....	103
6.6.4	Acquisition.....	105
6.6.5	Processing .....	105
6.7	1D Selective NOESY Experiment Using the Off- Resonance Option .....	107
6.7.1	Reference Spectrum .....	107
6.7.2	Selective Excitation Region Set Up .....	107
6.7.3	Setting Up the Selective NOESY .....	109
6.7.4	Acquisition.....	111
6.7.5	Processing .....	111
<b>7</b>	<b>2D Experiments using Shaped Pulses .....</b>	<b>113</b>
7.1	2D Edited HSQC Experiment with Adiabatic Pulses.....	113
7.1.1	Introduction .....	113
7.1.2	Sample.....	114
7.1.3	Reference spectrum.....	114
7.1.4	Setting up the HSQC experiment.....	114
7.1.4.1	Limit Setting .....	115
7.1.4.2	Acquisition.....	116
7.1.4.3	Processing .....	117
7.2	2D Selective HMBC Experiment.....	121
7.2.1	Introduction .....	121
7.2.2	Sample.....	121
7.2.3	Preparation Experiments .....	122
7.2.3.1	1D Proton Experiment.....	122
7.2.3.2	2D HMBC Experiment.....	122
7.2.3.3	1D Proton Decoupled Carbon Experiment.....	123
7.2.3.4	Method 1 for Setting Up the Selective HMBC Experiment.....	123
7.2.3.5	Selective Excitation Region Setup .....	124
7.2.3.6	Processing .....	127
7.2.3.7	Method 2 for Setting Up the Selective HMBC Experiment.....	127
7.2.3.8	Processing .....	130
<b>8</b>	<b>T1 Experiment.....</b>	<b>131</b>
8.1	Introduction .....	131
8.2	Sample.....	131
8.3	Proton Inversion-Recovery T1 Experiment .....	131
8.3.1	Preparation Experiment .....	131

8.3.2	Setting up the T1 Experiment .....	132
8.3.3	Limit Setting .....	133
8.3.4	Acquisition.....	135
8.3.5	Processing .....	135
8.3.6	T1 Calculation .....	136
<b>9</b>	<b>Pulse Calibration .....</b>	<b>143</b>
9.1	Introduction .....	143
9.2	Sample.....	143
9.3	$^1\text{H}$ 90° Transmitter Pulse.....	143
9.3.1	Preparation Experiment .....	143
9.3.2	Parameter Setup.....	144
9.3.3	Acquisition.....	146
9.3.4	Processing .....	146
9.3.5	Determine the $^1\text{H}$ 90° Transmitter Pulse .....	147
9.4	$^{13}\text{C}$ 90° Decoupler Pulse .....	150
9.4.1	Preparation Experiment .....	150
9.4.2	Parameter Setup.....	150
9.4.3	Determine the $^{13}\text{C}$ 90° Decoupler Pulse.....	153
9.5	$^{15}\text{N}$ 90° Decoupler Pulse .....	155
9.5.1	Parameter Setup.....	155
9.5.1.1	Two Channel System.....	157
9.5.1.2	Three Channel System .....	158
9.5.2	Determine the $^{15}\text{N}$ 90° Decoupler Pulse.....	159
<b>10</b>	<b><math>^1\text{H}</math> Homonuclear Decoupling .....</b>	<b>161</b>
10.1	Sample.....	161
10.2	$^1\text{H}$ Homonuclear Decoupling Experiment.....	161
10.2.1	Introduction .....	161
10.2.2	Acquisition.....	161
10.2.3	Processing .....	164
10.2.4	Plotting Two Spectra on the Same Page .....	165
<b>11</b>	<b>Proton DOSY Experiment .....</b>	<b>167</b>
11.1	Introduction .....	167
11.2	Sample.....	167
11.3	DOSY Experiment.....	167
11.3.1	Pulse Programs .....	167
11.3.2	Preparation Experiment .....	168
11.3.3	Acquisition.....	170
11.3.4	Processing .....	170
11.3.5	Limit Settings .....	170
11.3.6	Parameter Setup.....	171
11.3.7	Running the Experiment .....	175
11.3.8	Processing .....	178
11.3.9	Calculating the Diffusion Coefficient .....	179
11.3.10	Displaying the DOSY Plot.....	185
<b>12</b>	<b>Multiplet Analysis .....</b>	<b>187</b>
12.1	Introduction .....	187

12.2	Sample.....	187
12.3	Multiplet Assignments .....	187
12.3.1	Preparation Experiment .....	187
12.3.2	Limit Settings .....	187
12.3.3	Acquisition.....	187
12.3.4	Processing .....	188
12.3.5	Peak Picking .....	190
12.3.6	Assigning the Multiplets .....	191
<b>13</b>	<b>Adding a New Solvent.....</b>	<b>201</b>
13.1	Introduction .....	201
13.2	Adding Cyclohexane-d12 to the Solvent List .....	201
13.3	TopShim Solvent Parameters .....	207
<b>14</b>	<b>Contact .....</b>	<b>209</b>





# 1 About This Manual

This manual enables safe and efficient handling of the device.

This manual is an integral part of the device, and must be kept in close proximity to the device where it is permanently accessible to personnel. In addition, instructions concerning labor protection laws, operator regulations tools and supplies must be available and adhered to.

**Before starting any work, personnel must read the manual thoroughly and understand its contents.** Compliance with all specified safety and operating instructions, as well as local work safety regulations, are vital to ensure safe operation.

The figures shown in this manual are designed to be general and informative and may not represent the specific Bruker model, component or software/firmware version you are working with. Options and accessories may or may not be illustrated in each figure.

## 1.1 Symbols and Conventions

Safety instructions in this manual and labels of devices are marked with symbols.

The safety instructions are introduced using indicative words which express the extent of the hazard.

In order to avoid accidents, personal injury or damage to property, always observe safety instructions and proceed with care.

### DANGER



**DANGER:** Indicates a hazardous situation that, if not avoided, will result in death or serious injury. This signal word is limited to the most extreme situations.

This is the consequence of not following the warning.

1. This is the safety condition.
  - ▶ This is the safety instruction.

### WARNING



**WARNING:** Indicates a hazardous situation that, if not avoided, could result in death or serious injury.

This is the consequence of not following the warning.

1. This is the safety condition.
  - ▶ This is the safety instruction.



### CAUTION

**CAUTION:** Indicates a hazardous situation that, if not avoided, could result in minor or moderate injury.

This is the consequence of not following the warning.

1. This is the safety condition.
  - ▶ This is the safety instruction.

### NOTICE

**NOTICE:** Indicates information considered important, but not hazard-related (e.g. messages relating to property damage).

This is the consequence of not following the notice.

1. This is a safety condition.
  - ▶ This is a safety instruction.

### SAFETY INSTRUCTIONS

**SAFETY INSTRUCTIONS** are used for control flow and shutdowns in the event of an error or emergency.

This is the consequence of not following the safety instructions.

1. This is a safety condition.
  - ▶ This is a safety instruction.



This symbol highlights useful tips and recommendations as well as information designed to ensure efficient and smooth operation.

## 2 Introduction

### 2.1 General

---

This manual was written for AVANCE systems running **TopSpin version 4.x including patches** and should be used as a guide through the set up process for some experiments. The success of running the experiments in this manual is under the assumption that all parameters have been entered into the prosol table.

This manual features various advanced procedures for  $^1\text{H}$ ,  $^{13}\text{C}$ ,  $^{15}\text{N}$  and  $^{31}\text{P}$  experiments. It is assumed that the user is already familiar with acquisition and processing of simple 1D NMR spectra, chapter *1D Proton Experiment* and chapter *1D Carbon Experiments* described in the TopSpin Guide Book *Basic NMR Experiments*.

### 2.2 Disclaimer

---

This guide should only be used for its intended purpose as described in this manual. Use of the manual for any purpose other than that for which it is intended is taken only at the users own risk and invalidates any and all manufacturer warranties.

Some parameter values, especially power levels suggested in this manual may not be suitable for all systems (e.g. Cryo probes) and could cause damage to the unit. Therefore only persons trained in the operation of the AVANCE systems should operate the unit.

#### **NOTICE**


##### **Material Damage Due to Excessive Power**

The NMR probe can be severely damaged if too much power or power over a too long time is applied.

- ▶ Always start to optimize pulses with low power values and short pulses. Respect the pulse and power limits as programmed into the PICS data of the probe.



## 3 The TopSpin Interface



1	Menu Bar	8	Status Display Bar
2	Workflow Button Bar	9	Dataset Window
3	Tool Bar	10	Dataset Window Tabs
4	Browser and Search Window	11	Print, Export, Copy and Publish
5	Structure Window, Command Line History, Status Line History	12	Viewing Options
6	Command Line	13	Window Switcher, Login, Setup Preferences and Help
7	Current Dataset Bar		


### Setup Preferences

For all changes in the TopSpin appearance use the global **Setup preferences** button in the menu bar.

TopSpin can be tailored to your preference in many respects. This ranges from startup options to spectrum objects, menu settings, remote connections, colors and fonts etc. Every standard user can create his own set of preferences.

A dialog box will appear with, at the left side, the categories that can be tailored. Click the category you want to view or change. It will become highlighted and the corresponding objects will be displayed at the right part of the dialog box.


# The TopSpin Interface



## TopSpin on High Resolution Screens


To adapt the font and icon size to a small/standard/medium or large screen resolution

- click **Setup Preferences** and in the category **Window Settings | All In One Forts & Icon Size** click **Change**.




## Help/About TopSpin/Version and License Information

This button gives information about the TopSpin documentation, software version and license.




Clicking **Manuals (docs)** or entering **docs** on the command line will open the list of all manuals delivered with TopSpin.

- Click **Help | Manuals | Acquisition Application Manuals | Dosy** will open the DOSY manual, for example.




# The TopSpin Interface

Enter **help** to get information for an individual command. Three different sources can be selected for the search:



## Tooltips

Pointing to a button with the mouse in the various menus opens a tooltip that describes the button functionality. Example:




## Workflow Button Bar

The workflow-based interface with its arrangement of all working processes allows the user to control the workflow intuitively.

Clicking one of the menu buttons opens the corresponding workflow. It contains a horizontal feature list which stays open and provides all functionality for this workflow with one mouse-click.


Furthermore, some of the buttons on the Workflow button bar include a **drop-down** arrow. Click the **drop-down** arrow to see more options.






## Dataset Toolbar

Depending on which dataset window tab is selected, an individual dataset toolbar is displayed, in the example the ACQU toolbar:



**Note:** The ACQU window tab is only displayed when TopSpin controls a spectrometer (noticeable through the Acquire tab in the TopSpin menu). When TopSpin is installed for processing-only, the Acquire tab is also not displayed.


**Note:** Sometimes the dataset toolbar is displayed without dataset window tabs, e.g. entering a command as `.bas1` will display the baseline correction toolbar:



## Browse and Search Window


The Browser window provides tabs as:

- Data browser and Search
- History browser
- Dataset Switcher
- Experiment Selector library



# The TopSpin Interface


The window can be toggled On or Off with a click on the black **left** or **right** arrow.

Alternatively, the **Toggle Perspective** button on the top-right of the TopSpin window can be clicked. 

The data tree tab includes a search field. Enter a search phrase, for example:


Search:  X

To limit the search results, select or deselect the data trees with **SHIFT** click or **CTRL** click.





## History Tab

The History tab displays the last opened datasets in a list.




## Dataset Switcher Tab

This tab has a similar function as the Window Switcher , see below, but only provides a quick overview of all currently opened datasets.




## Experiment Library Tab

The Experiment library tab provides easy access to the vast experiment library of standard experiments that is available in the standard release of TopSpin. It is also a tool that allows the user to personalize his most frequently used experiments into a separate library.



## Structure Window

Molecular structures as *.mol* or *.pdb* files can be displayed here and are freely resizable. The structure window can be toggled On or Off with a click on the black **up** or **down** arrow ▲▼.



- Drag the vertical or horizontal split bar to resize the structure window.


## Window Switcher



The **Window Switcher** button indicates the number of all opened dataset windows plus all opened TopSpin window types like

- Dataset windows.
  - Lock display window.
  - Acquisition display window.
  - BSMS display window.
  - Temperature unit window.
- Click the **Window Switcher** button to switch between these windows.



A dataset can be closed with the **Close** button .

## Viewing Integrals, Peaks and other Spectra Components


- Click **Spectrum display properties** to toggle the visibility of integrals, peaks and other spectra components. They can only be displayed when available.



## Shift, Scroll and Zoom Spectra Axes

Clicking and dragging the spectrum axis allows intuitive scrolling through the spectrum.

- Single Clicks on the end regions of the axes or a click and dragging the mouse shift through the axes depending on the mouse position, see the next figure.
- Double click in the middle of the spectra axes switches to full axes region.
- Double click in the middle of the spectra itself maximizes all axes regions.
- Click on axes to shift left, right, up or down.
- Turn mouse wheel to zoom in or out.



## Dataset Windows


The TopSpin window has a dataset area that may contain multiple dataset windows.


To open multiple dataset windows

- Right-click in the Browser window and in the list, select **Display in new window**.

Note, that selecting **Display** will override a current dataset.

The size of the data area depends on the overall size of the TopSpin window and on presence of the Browser. Note, that the Browser window can be toggled On or Off with the **Tog-**

**gle Perspective** button . The following figure shows the TopSpin window without the Browser and three dataset windows.




# The TopSpin Interface

## How to Arrange Dataset Windows

If the data area contains multiple dataset windows, you can arrange them in various ways. All the arrange commands arrange the windows left to right and/or top to bottom in the order in which the windows have been active. The currently active dataset window will therefore be positioned at the top and/or left of the data area.


To arrange the dataset windows as a grid:

- Click **Show Layout Options**  and **Show as Grid** .

Depending on the number of windows, they will be arranged vertically and/or horizontally.


To arrange dataset windows in stack:

- Click **Show Layout Options**  and .



To arrange dataset windows side by side:

- Click **Show Layout Options**  and .



To display a dataset windows as full screen:

- Click **Show Layout Options**  and  or click the **full screen** windows button .

To close the active dataset window:

- Click **File | Close Active Window** or enter **Crtl-w**.  
or



- Click the  button in the windows title bar.

To close all dataset windows:

- Click **File | Close All Windows** or enter **closeall**.


## How to Swap Dataset Windows

Within a certain layout, you can easily swap two TopSpin windows with the command **swin**. If the data area contains exactly two windows, **swin** simple swaps their position and geometry.

If it contains more than two dataset windows, **swin** opens a list from which you can select any window to be swapped with the currently selected (active) window.


## Acquisition Status Bar

- Right-click the Status display or Current dataset bar to toggle the Acquisition status bar On or Off.




- Click **Status Bar Preferences** to set the *Auto-open the acquisition status bar* option as default.


Acquisition Status Bar On/Off  
Status Bar Preferences



The acquisition status bar contains the new Spectrometer Status area to turn the spectrometer On and Off. This functionality is also available in the menu with a click on **Manage | Spectrometer | Spectrometer power On/Off** or the command **pdudisp**.



## Print, Export and Publish

Click on any of the 4 icons  to either print the active window, export the active window to a PDF file, copy the active window to the clipboard or, show more publishing options such as E-mailing a dataset, sending a dataset to cloud storage, etc.

- Click **Show More publishing options**, e.g. to copy and paste, E-mail or use shared cloud directories.




## 4 1D Solvent Suppression Experiments

### 4.1 Introduction

---

Many experiments on samples dissolved in protonated solution require some method to minimize the strong resonance belonging to the solvent. This suppression can be performed in several ways, depending on the number of signals to suppress, and on which part of the pulse sequence can be modified. Solvent suppression can be applied during the relaxation period just prior to the conventional pulse sequence as outlined in the figure below. This is referred to as Pre-saturation.



However, pre-saturation can also reduce the signal intensities of exchangeable protons. For this reason, other schemes, as the WATERGATE, WET and Excitation Sculpting schemes, can be used to overcome this problem and are discussed in this chapter.

In HPLC-NMR applications it is mandatory to suppress multiple-solvent resonances. The incorporation of specific multiple-solvent suppression schemes into pulse sequences is made in analogy with classical methods.

### 4.2 Samples

---

2 mM Raffinose in 90% H<sub>2</sub>O + 10% D<sub>2</sub>O

2 mM Lysozyme in 90% H<sub>2</sub>O + 10% D<sub>2</sub>O

### 4.3 Preparation Experiment

---

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the Create New Dataset window, enter or select:
  - NAME = **solvent\_suppression\_exp**
  - EXPNO = 1
  - Directory = e.g. C:\DMB



For the following steps, use the Workflow button bar from left to right.

- Click **Sample** and eject the sample, if there is one inserted, and insert the new sample.
- Click **Lock** and select **H2O + D2O** solvent.
- To tune the probe, click **Tune**.
- Click **Spin** and select **Turn sample rotation off**.



Solvent suppression experiments should be run non-spinning.

- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent depended parameters, click **Prosol**.

## 4.3.1 Acquisition


- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.


## 4.3.2 Processing

- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.

This executes a processing program including commands such as an exponential window function **em**, Fourier transformation **ft**, an automatic phase correction **apk** and a baseline correction **abs**.


- On the **Proc. Spectrum** button, click the drop-down arrow to see more options.
- In the list, select **Configure Standard Processing (proc1d)**.







## 4.4 1D Solvent Suppression with Pre-saturation

**Pre-saturation** is the most common procedure to minimize and suppress the intense solvent resonance when <sup>1</sup>H spectra are recorded in protonated solutions. This experiment is performed by applying a low-power continuous wave irradiation on the selected resonance during the pre-scan delay:




### 4.4.1 Parameter Setup

- On the command line, type:  
**wrpa 2**  
**re 2**
- Expand the peak at **4.7ppm**.
- On the toolbar, click **Set RF from cursor**. 



- Move the cursor line to the center of the peak and press the mouse button.




- In the O1/O2/O3 window, click **O1**.
- In the Dataset window, select the **Spectrum** tab.


## 4.4.2 Fine Tuning

- On the menu bar, click **Acquire**.
- On the Workflow button bar, click **Gain**.
- On the **Run** button, click the **drop-down** arrow to see more options.
- In the list, select **Real-Time Go Setup (gs)**.

Start Acquisition (zg)
Start acquisition, add to existing data (go)
Transfer FID To Disk (tr)
Estimate Exp. Time (expt)
Real-Time Go Setup (gs)
Optimize Acquisition Params (popt)
Start Automation AU program (xaau)

# 1D Solvent Suppression Experiments

- Click **Unshuffle data on display**. 
- In the Go Setup (gs) window, select the **Offset** tab.




- Change the **O1** value by clicking just below or above the adjust slider.




For smaller changes, adjust the **sensitivity** inside the gs window to a smaller value.

- Optimize the fid area in the Acquisition information window for a smaller integration value and for the FID to become a single exponential decay.
- In the Go Setup (gs) window, click **Save**.
- In the Go Setup (gs) window, click **Stop**.





- In the Save changed GS parameters window, click **OK**.




## 4.4.3 Acquisition

- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.


## 4.4.4 Processing

- Process and phase correct the spectrum.

## 1D Solvent Suppression Experiments



The figure above shows the solvent suppressed 1D spectrum of the Raffinose sample. The figure below shows the 1D spectrum of the Lysozyme sample.



## 4.5 1D Solvent Suppression with Presaturation and Composite Pulses

---

This experiment is performed by applying a low-power continuous wave irradiation on the water resonance during the pre-scan period, followed by a rapid succession of four 90° pulses to further reduce the residual hump of the water signal:



### 4.5.1 Parameter Setup

---

- Follow the instructions *Parameter Setup* and *Fine Tuning* in the chapter [1D Solvent Suppression with Pre-saturation \[▶ 30\]](#).
- In the Dataset window, select the **AcquPars** tab.
- Enter:  
PULPROG = **zgcppr**



There is also a standard parameter set **ZGCPPR**.

---

- In the Dataset window, select the **Spectrum** tab.

### 4.5.2 Acquisition

---


- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

### 4.5.3 Processing


---

- Process and phase correct the spectrum.

## 1D Solvent Suppression Experiments




The figure above shows the solvent suppressed 1D spectrum of the Raffinose sample. The following figure shows the 1D spectrum of the Lysozyme sample:




### 4.6 1D Solvent Suppression Using the Noesy Sequence

This experiment is performed by using the 1D version of the noesyphpr sequence applying a low-power continuous wave irradiation on the water resonance during the pre-scan and during the mixing time period of the NOESY sequence:



### 4.6.1 Parameter Setup

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- Enter or select  
 Name: **noesypr1d**  
 EXPNO: 1  
 Read parameterset: **NOESYPR1D**  
 Set solvent: **H2O+D2O**




- On the Workflow button bar, click **Lock** and select **H2O+D2O**.
- On the Workflow button bar, click **Tune**.
- On the Workflow button bar, click **Shim**.
- On the Workflow button bar, click **Gain**.
- On the Workflow button bar, click **Run**.


### 4.6.2 Processing

- Process and phase correct the spectrum.

## 1D Solvent Suppression Experiments




The figure above shows the solvent suppressed 1D spectrum of the Raffinose sample. The figure below shows the 1D spectrum of the Lysozyme sample.



## 4.7 1D Solvent Suppression with WATERGATE


The WATERGATE (**WATER** suppression by **Gr**Adient **Tai**lored **Exc**itation) technique, which uses pulsed field gradients, is claimed to be independent of line-shape, yielding better suppression compared with other methods. Exchangeable protons are not affected and there is no phase jump at the water resonance, although signals very close to the water resonance are also suppressed.

The sequence is in principle, a spin-echo experiment in which the 180 degree pulse is embedded between two pulsed field gradients. After excitation by the first pulse p1 the field gradient G1 dephases all coherence. The selective inversion element consists of a symmetrical 3-9-19 pulse sequence 3a-t-9a-t-19a-t-19a-t-9a-t-3a, with  $2\theta = 180^\circ$ , as shown in the figure below. Additional suppression appears at different sidebands (1/t).



### 4.7.1 Parameter Setup


- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- Enter or select  
 Name: **p3919gp**  
 EXPNO: 1  
 Read parameterset: **P3919GP**  
 Set solvent: **H2O+D2O**



- On the Workflow button bar, click **Lock** and select **H2O+D2O**.
- On the Workflow button bar, click **Tune**.
- On the Workflow button bar, click **Shim**.
- On the Workflow button bar, click **Gain**.
- On the Workflow button bar, click **Run**.

### 4.7.2 Processing

- Process and phase correct the spectrum.




The figure above shows the solvent suppressed 1D spectrum of the Raffinose sample. The figure below shows the 1D spectrum of the Lysozyme sample.




The final pulse, p0, of the 3-9-19 pulse train can be optimized (shortened) in gs mode to eliminate the dispersive tail on the HOD peak.






## 4.8 1D Solvent Suppression with Excitation Sculpting



### 4.8.1 Parameter Setup

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- Enter or select  
 Name: **zgesgp**  
 EXPNO: 1  
 Read parameterset: **ZGESGP**  
 Set solvent: **H2O+D2O**


# 1D Solvent Suppression Experiments




- On the Workflow button bar, click **Lock** and select **H2O+D2O**.
- On the Workflow button bar, click **Tune**.
- On the Workflow button bar, click **Shim**.
- On the Workflow button bar, click **Gain**.
- On the Workflow button bar, click **Run**.

## 4.8.2 Processing

- Process and phase correct the spectrum.




The figure above shows the solvent suppressed 1D spectrum of the Raffinose sample. The figure below shows the 1D spectrum of the Lysozyme sample.



## 4.9 1D Solvent Suppression with WET

This pulse sequence uses a shaped, selective pulse and pulse field gradients to suppress one or more solvent signals. The option of carbon decoupling is available for suppression of solvent signals with large  $^{13}\text{C}$  satellites. It provides very efficient suppression with excellent selectivity.



### 4.9.1 Sample

2 mg Sucrose in Acetonitrile and D<sub>2</sub>O

## 4.9.2 Preparation Experiment

---

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the Create New Dataset window, enter or select:  
NAME = **wet\_solvent\_suppression\_exp**  
EXPNO = **1**  
Experiment: Select **PROTON**  
Set Solvent: Select **CH3CN+D2O**
- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.

For the following steps, use the Workflow button bar.

- Click **Sample** and **Insert**.
- Click **Lock** and select the **CH3CN+D2O** solvent.
- Click **Tune** to tune the probe.
- Click **Spin** and select **Sample rotation off**.



Solvent suppression experiments should be run non-spinning.

---

- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent depended parameters, click **Prosol**.

## 4.9.3 Acquisition

---

- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.


## 4.9.4 Processing

---

- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.


This executes a processing program including commands such as an exponential window function **em**, Fourier transformation **ft**, an automatic phase correction **apk** and a baseline correction **abs**.

- On the **Proc. Spectrum** button, click the drop-down arrow to see more options.
- In the list, select **Configure Standard Processing (proc1d)**.



## 4.9.5 Selective Excitation Region Setup

- On the menu bar, click **Acquire**.
- On the **More** button, click the **drop-down** arrow to see more options.




- In the list, select **Setup Selective 1D Expts.**

The Workflow button bar changes for setting up the 1D selective experiment.

- On the Workflow button bar, click **1D Selective Experiment Setup**.




# 1D Solvent Suppression Experiments



There is no other function to this button then the instruction displayed above.


- In the message window, click **Close**.
- On the Workflow button bar, click **Define Regions**.




The Define Regions toolbar is displayed:




- Integrate the regions on the peaks to be suppressed (e.g. **4.7 ppm** and **2.45 ppm**).




- Click **Save / Export integration regions.** 
- Select **Save Regions to 'reg'**.



- To exit the integration mode, click **Return do NOT save regions!** 
- In the Save Changes window, click **No**.




- On the **Create Dataset** button, click the **drop-down** arrow to see more options.
- In the list, select **Mult. Solvent Suppr./WET**.




- In the LC1DWTDC window, enter NS = **16**
- In the LC1DWTDC window, click **Accept**.

# 1D Solvent Suppression Experiments




- Check the parameters in the information window.



- In the sel1d window, click **OK** to start the acquisition.


## 4.9.6 Processing

- On the menu bar, click **Process**.
- On the **Proc Spectrum** button, click the **drop-down** arrow to see more options.
- Select **Configure Standard Processing (proc1d)**.




- In the proc1d window deselect the following options:
  - **Auto-Phasing (apk)**
  - **Set Spectrum Reference (sref)**
  - **Auto-Baseline correction (abs)**
  - **Warn if Processed data exist**





- In the proc1d window, click **Execute**.




- On the Workflow button bar, click **Adjust Phase**.

The Dataset window tabs are replaced with the **Adjust Phase** toolbar.




- Move the cursor line on top of the peak at **5.7 ppm** and right-click.
- On the shortcut menu, select **Set Pivot Point**.




- Adjust the **0**-order phase on the selected pivot point to positive absorption.
- Adjust the **1<sup>st</sup>** order phase on the peaks between **4.6 ppm** and **3.6 ppm**.

# 1D Solvent Suppression Experiments



The peak at **5 ppm** and **2.4 ppm** are the suppressed solvent peaks and those will appear out of phase.

- To store the phase values, click **Return and Save phased spectrum.**



# 5 2D Gradient Experiments

## 5.1 Introduction

The vital importance of NMR in chemistry and biochemistry relies on the direct relationship between any given NMR experiment and the molecular information that can be extracted from it. Thus, every experiment is based on some NMR parameter, usually coupling constants or NOE, which is related to a specific molecular parameter (through-bond or through-space connectivity, chemical exchange, molecular motion...). The quantitative measurement of such NMR parameters allows us to obtain valuable information about structural parameters such as dihedral angles, intermolecular distances, relaxation and exchange rates. etc... For this reason, the development of new and/or improved NMR methodologies is a key factor to be considered. Since the 1990's when the gradients were introduced as a useful tool to incorporate them into NMR applications, the suite of NMR experiments available to researchers has grown. A large percentage of them are using pulse field gradients.

Gradient enhanced NMR spectroscopy is widely used in liquid state spectroscopy for coherence pathway selection, solvent suppression, artifact reduction, and diffusion weighting and has had a tremendous impact by improving the quality of NMR spectra.

Thus, all advantages offering the incorporation of PFG (Pulsed Field Gradients) as a powerful element into high-resolution NMR pulse sequences, combined with the advanced software tools available at the present time to acquire and process multidimensional NMR experiments with great simplicity, has dramatically changed the concept of routine work in NMR for chemists.


## 5.2 Setting up the COSY Experiment



This section illustrates the need for phase cycling to eliminate spectral artifacts in a non-gradient experiment.

The steps below assume that the sample remains in the magnet after observing the Proton spectrum.

The 2D COSY (**C**ORrelated **S**pectroscop**Y**) is a two dimensional homonuclear experiment that correlates through bond coupled protons by tracing out connectivity via the homonuclear  $J_{HH}$  coupling constant. The COSY experiment can be acquired with or without PFG. The PFG selects coherences for observation and eliminates spectral artifacts, whereas in the non-gradient version, this is accomplished by phase cycling of the RF pulses in the sequence.



- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
  - In the New Dataset window, enter: NAME = **cosyqf\_exp**
  - EXPNO = **1**
  - Experiment: select **COSY90SW**
  - Set Solvent: select **DMSO**
  - Set Title: **non-gradient, magnitude mode COSY, NS=1**
  - In the New Dataset window, click **OK**.
  - Set TD(F1) to **128**.
  - Set NS to **1**.
- Under the Acquire menu bar, use the Workflow button bar.
  - Click **Tune** to tune the probe.
  - Click **Spin** and select **Sample rotation off**.



---

2D experiments should be run non-spinning.

---

- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent dependent parameters, click **Prosol**.

### 5.2.1 Limit Setting

---


- On the Workflow button bar, click **SetLimits**.
- To open the 1D Proton spectrum, right click on the dataset name in the browser window (e.g. proton\_exp) and select **Display** or click and hold the left mouse button for dragging the 1D Proton dataset into the spectrum window.
- Expand the spectrum to display all peaks, leaving ca. **1.0 ppm** of baseline on either side of the spectrum.




---

The solvent peak may be excluded if it falls outside of the region of interest. Digital filtering however is only applied in F2 and the solvent peak will be folding in F1.

---



- In the SetLimits message window, click **OK** to assign the new limit.
- In the message window click **Close**.



The display changes back to the 2D dataset.

## 5.2.2 Acquisition


- To auto-adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

## 5.2.3 Processing

- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.


## 2D Gradient Experiments

This executes a standard processing program **proc2d**. To configure this program or select the right options, click the down arrow inside the **Proc. Spectrum** button. Since this is a magnitude mode experiment the phase correction **apk2d** should be disabled.




To eliminate the anti-diagonal and axial peaks using the phase cycling coded in the pulse program, the number of scans needs to be increased.

- Type the command **ixpno** on the command line to copy the parameters of the current dataset to a new dataset whose EXPNO is increased by 1.
- Set NS to **4**.
- Set Title: non-gradient, magnitude mode COSY, NS=4
- On the menu bar, click **Acquire | Run**.
- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.



### 5.3 Setting up the COSYGP Experiment

This section illustrates the recording of N- or P-type coherence depending upon the gradient values.



The steps below assume that the sample remains in the magnet after observing the proton spectrum.

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the New Dataset window, enter or select: NAME = **cosygp\_exp**
  - EXPNO = 1
  - Experiment: select **COSYGPSW**
  - Set Solvent: select **DMSO**
  - Set Title: gradient, magnitude mode COSY, N-type
- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.

For the following steps, use the Workflow button bar.

- Click **Tune** to tune the probe.
- Click the Spin down arrow and select **Sample rotation off**.




2D experiments should be run non-spinning.

- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent depended parameters, click **Prosol**.


### 5.3.1 Limit Setting

- On the Workflow button bar, click **SetLimits**.
- To open the 1D Proton spectrum, right click on the dataset name in the browser window (e.g. proton\_exp) and select **Display** or click and hold the left mouse button for dragging the 1D Proton dataset into the spectrum window.
- Expand the spectrum to display all peaks, leaving ca. **1.0 ppm** of baseline on either side of the spectrum.

The solvent peak may be excluded if it falls outside of the region of interest. Digital filtering however is only applied in F2 and the solvent peak will be folding in F1.



- In the SetLimits message window, click **OK** to assign the new limit.
- In the message window, click **Close**.



The display changes back to the 2D dataset.




### 5.3.2 Acquisition

- To auto-adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

### 5.3.3 Processing

- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.

This executes a standard processing program **proc2d**. To configure this program or select the right options, click the down arrow inside the **Proc. Spectrum** button. Since this is a magnitude mode experiment the phase correction **apk2d** should be disabled.




### 5.3.4 N- to P-Type Coherence Selection

To select the P-type coherence pathway instead of the N-type, the two Z gradient pulses in the pulse program need to have opposite sign.

- Type the command **ixpno** on the command line to copy the parameters of the current dataset to a new dataset whose EXPNO is increased by 1.

- Edit the pulse program **cosygpppqf** by changing the second occurrence of p16:gp1 in the sequence to p16:gp1\*-1. Save the edited pulse program to **cosygpppqf.Ptype** in the user subdirectory of the standard pulse program library.
- Change the PULPROG to **cosygpppqf.Ptype**.
- Set Title: gradient, magnitude mode COSY, P-type
- On the menu bar, click **Acquire | Run**
- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.




### 5.4 Setting up the Multiple-Quantum Filtered COSY Experiment

The COSY Multiple-Quantum Filtered (COSY-MQF) experiment is an alternative version of the COSY experiment, in which a multiple-quantum filter is inserted to allow the detection of signals from all coupled spin systems but suppresses signals arising from lower coherence levels. Thus, a COSY with a double-quantum filter (2D COSY-DQF experiment) efficiently suppresses single-quantum coherence from singlet uncoupled signals as, for instance, those of methyl groups or solvents. The COSY-DQF experiment can be performed in magnitude or phase-sensitive mode by selecting the appropriate phase programs and transform algorithm. However, phase-sensitive data is usually recommended.

In spectrometers equipped with gradient technology, gradient-based COSY versions are highly recommended.


The gp-2D COSY-MQF experiment yields a 2D COSY-MQF spectrum with a single scan per t1 increment provided that the S/N ratio is adequate. The main advantage of such an approach is the large reduction in the total acquisition time compared with a conventional phase-cycled 2D COSY-MFQ experiment. Magnitude-mode (**cosygpmpfqf**) or phase-sensi-

tive (**cosygpmfph**) data is obtained depending on the selected pulse sequence and acquisition/processing procedure. The COSY-MQF experiment traces out through-bond proton-proton connectivity via the homonuclear  $J_{\text{HH}}$  coupling constant.



### 5.4.1 Preparation Experiment

Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment*, *Experiment Setup* through *Processing*.



### 5.4.2 Setting up the MQF-COSY Experiment

The steps below assume that the sample remains in the magnet after observing the proton spectrum.

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the New Dataset window, enter or select:
  - NAME = **cosydqf\_exp**
  - EXPNO = **1**
  - Experiment: select **COSYGPDPHPSW**
  - Set Solvent: select **DMSO**
- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.

For the following steps, use the Workflow button bar.

- Click **Tune** to tune the probe.
- Click **Spin** and select **Sample rotation off**.




2D experiments should be run non-spinning.

- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent depended parameters, click **Prosol**.


### 5.4.3 Limit Setting

- On the Workflow button bar, click **SetLimits**.
- To open the 1D Proton spectrum, right click on the dataset name in the browser window (e.g. proton\_exp) and select **Display** or click and hold the left mouse button for dragging the 1D Proton dataset into the spectrum window.




- Expand the spectrum to display all peaks, leaving ca. **1.0 ppm** of baseline on either side of the spectrum.

The solvent peak may be excluded if it falls outside of the region of interest. Digital filtering however is only applied in F2 and the solvent peak will be folding in F1.



- In the setlimits message window, click **OK** to assign the new limit.
- In the message window click **Close**.



The display changes back to the 2D dataset.

#### 5.4.4 Acquisition


The first increment of the DQF-COSY experiment has a low signals to noise ratio and the signals grow as the experiment is progressing. It is therefore not advisable to use the automatic receiver gain adjustment **rga** since it adjusts the receiver gain on the first increment. In this case an AU program **au\_zgcosy** is available. Executing this AU program changes the pulse program to **zg** and performs an **rga** and then changes back again to **cosygpmfph** and then starts the acquisition.

- At the command prompt, type **au\_zgcosy**.


### 5.4.5 Processing


- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.

This executes a standard processing program **proc2d**. To configure this program or select the right options, click the down arrow inside the **Proc. Spectrum** button. Since this is a phase sensitive experiment the phase correction **apk2d** should be enabled.




### 5.4.6 Plotting

Click on any of the 4 icons  to either print the active window, export the active window to a PDF file, copy the active window to the clipboard or, show more publishing options such as E-mailing a dataset, sending a dataset to cloud storage, etc.



If desired, any changes can be administered by using the tools on the left side of the display.

- In the Print section, click the **down** arrow button and select **Print**.




## 5.5 Setting up the $^{13}\text{C}$ -HMBCGP Experiment

The **2D gradient HMBC** experiment records qualitative heteronuclear long-range connectivity, including through hetero nuclei. This section of the manual will guide you through the setup of an  $^1\text{H}/\text{X}$  gradient experiment using the standard Bruker HMBCGP parameter set. In

## 2D Gradient Experiments

In addition to changing the nucleus in F1 from  $^{13}\text{C}$  to another X-nucleus, the gradient ratio for the new X-nucleus also must be calculated. The HMBC pulse sequence is shown in the figure below.



The time intervals depicted in the pulse sequence diagrams are not drawn to scale. For example,  $d1$  is typically a few seconds while  $p1$  is typically a few microseconds in length.

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the New Dataset window, enter or select: NAME = **13Chmbcgp\_exp**
  - EXPNO = 1
  - Experiment: select **HMBCGP**
  - Set Solvent: select **DMSO**
  - Set Title: gradient, magnitude mode HMBC
- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.

For the following steps, use the Workflow button bar.

- Click **Tune** to tune the probe.
- Click the Spin down arrow and select **Sample rotation off**.



2D experiments should be run non-spinning.


- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent depended parameters, click **Prosol**.




### 5.5.1 Limit Setting

- On the Workflow button bar, click **SetLimits**.
- To open the 1D Proton spectrum, right click on the dataset name in the browser window (e.g. proton\_exp) and select **Display** or click and hold the left mouse button for dragging the 1D Proton dataset into the spectrum window.
- Expand the spectrum to display all peaks, leaving ca. 1.0 ppm of baseline on either side of the spectrum.

The solvent peak may be excluded if it falls outside of the region of interest. Digital filtering however is only applied in F2 and the solvent peak will be folding in F1.



- In the SetLimits message window, click **OK** to assign the new limit.
- In the message window click **Close**.



The display changes back to the 2D dataset.


### 5.5.2 Acquisition

- To auto-adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

### 5.5.3 Processing

- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.

This executes a standard processing program **proc2d**. To configure this program or select the right options, click the down arrow inside the **Proc. Spectrum** button. Since this is a magnitude mode experiment the phase correction **apk2d** should be disabled.




## 5.6 Setting up the 15N-HMBCGP from the 13C-HMBCGP

Although there is a dedicated standard parameter set for acquiring a 15N-HMBC, HMBCGP\_15N, in this section one will be created starting from the previously acquired 13C-HMBC. In the new dataset, the X nuclei will be changed from 13C to 15N. This exercise illustrates that in the gradient HMBC experiment the gradients are used to select the coherence corresponding to the HX nuclei defined in the dataset. These steps can be used to set up the HMBC for other heteronuclei combinations, e.g. 1H/31P, for which there is no dedicated parameter set in the library.

### 5.6.1 Limit Setting


- Type the command **ixpno** on the command line to copy the parameters of the current dataset to a new dataset whose EXPNO is increased by 1.

- Change the Title to **15Nhmbcgp\_exp**
- On the Workflow button bar, click **SetLimits**.
- To open the 1D Proton spectrum, right click on the dataset name in the browser window (e.g. proton\_exp) and select **Display** or click and hold the left mouse button for dragging the 1D Proton dataset into the spectrum window.




- Expand the spectrum to display all peaks, leaving ca. **1.0 ppm** of baseline on either side of the spectrum.


The solvent peak may be excluded if it falls outside of the region of interest. Digital filtering however is only applied in F2 and the solvent peak will be folding in F1.



- In the setlimits message window, click **OK** to assign the new limit.
- In the message window click **Close**.

## 5.7 Setting up the 15N-HMBC Experiment

- In the Dataset window, select the AcqPars tab.
- Click the **Set nuclei and routing** button  to display the routing window.
- Change the following parameter using the pull-down arrow:  
NUC2 = **15N**





After switching NUC2 from 13C to 15N, click **Default** and then **Save and Close**.




Do not modify the routing!

After closing the Channel Routing window, two windows will appear to alert the user to update some experimental parameters in the AcqPars tab of the dataset:






- In the AcqPars tab scroll down to the Nucleus 2 section. Either ensure that 15N is defined as the 2<sup>nd</sup> nucleus or select it using the Edit button.
- Set O2P [ppm] to **150**.



- Click on the getprosol icon at the top of the AcqPars tab to update the 15N 90 degree pulse:
- Change the following parameters:  
**PULPROG = hmbcgpndqf**  
**D1 = 2 sec**  
**CNST13 = 5 Hz**  
**SW [ppm] in F1 = 400 ppm**
- To calculate the gradient ratio necessary for 1H-15N detection, type **gradratio** on the command line. The correct gradient ratios will appear in a pop-up window and will be updated in the dataset automatically.



- Under the Acquire menu, click on **Tune** to tune and match the probe to 15N.

### 5.7.1 Acquisition


- To auto-adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

### 5.7.2 Processing


- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.

This executes a standard processing program **proc2d**. To configure this program or select the right options, click the down arrow inside the **Proc. Spectrum** button. Since this is a magnitude mode experiment the phase correction **apk2d** should be disabled.

- To correct the F1 axis for <sup>15</sup>N, click on **Calib. Axis** under the Process menu.




- At the top left in the dataset window, click on **C** to Calibrate to center of spectrum. Set F1[ppm] to 150 and click **OK**.



## 5.8 2D $^1\text{H}$ / $^{31}\text{P}$ Gradient HMBC Experiment


The **2D gradient HMBC** experiment records qualitative heteronuclear long-range connectivity, including through hetero nuclei. This section of the manual will guide you through the set up of a  $^1\text{H}/\text{X}$  gradient experiment using the standard Bruker HMBCGP parameter set. In addition to changing the nucleus in F1 from  $^{13}\text{C}$  to another X-nucleus, the gradient ratio for the new X-nucleus also has to be calculated. The HMBC pulse sequence is shown in the figure below.



The time intervals depicted in the pulse sequence diagrams are not drawn to scale. For example, d1 is typically a few seconds while p1 is typically a few microseconds in length.

### 5.8.1 Sample


30mg Dimethyl methylphosphonate in  $\text{CDCl}_3$



This  $^{31}\text{P}$  nucleus in this sample does not have any direct proton attached. The long range coupling from the 3 methyl protons to  $^{31}\text{P}$  is **17 Hz**, where the other 6 methyl protons through the additional oxygen nuclei show a J-value of **11 Hz**. The J-values can be easily obtained from the proton spectrum (see chapter  [\$^1\text{H}\$  Reference Experiment \[ 71 \]](#)).

### 5.8.2 $^1\text{H}$ Reference Experiment

Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book Basic NMR experiments, chapter *1D Proton Experiment*, Paragraph *Experiment Setup* through *Processing* using  $\text{CDCl}_3$  as a lock solvent.



### 5.8.3 <sup>31</sup>P Reference Experiment

The steps below assume that the sample remains in the magnet after observing the proton spectrum.

- On the menu bar, click **Acquire** | **Create Dataset** to open the Create New Dataset window.
- In the New Dataset window, enter or select:  
NAME = **31P\_exp**  
EXPNO = **1**  
Experiment: select **P31CPD**  
Set Solvent: select **CDCl3**
- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.

For the following steps, use the Workflow button bar.

- Click **Tune** to tune the probe.
- Click **Spin** and select **Sample rotation on**.
- Click **Shim** - for best homogeneity use TopShim.
- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- Click **Prosol** to load the probe/solvent depended parameters.




### 5.8.4 Acquisition


- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

### 5.8.5 Processing


- Process and phase correct the spectrum.



### 5.8.6 Setting up the HMBC Experiment

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the New Dataset window, enter or select:  
NAME = **1H\_31P\_hmbc\_exp**  
EXPNO = 1  
Experiment: **HMBCGP**  
Set Solvent: **CDCl3**
- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.
- In the Dataset window, select the **AcquPars** tab.
- Click the **Set nuclei and routing** button to display the routing window. 

## 2D Gradient Experiments



- Change the following parameter:  
NUC2 = 31P



Do not modify the routing!

- Click the **Save and Close** button inside the Channel Routing window.
- Scroll down to the Nucleus 1 section in the AcqPars window.

Nucleus 1				
NUC1	1H	Edit...	31P	Observe nucleus
O1 [Hz]	1200.68		5570.91	Transmitter frequency offset
O1P [ppm]	4.001		45.853	Transmitter frequency offset
SFO1 [MHz]	300.1312007		121.5004219	Transmitter frequency
BF1 [MHz]	300.1300000		121.4948510	Basic transmitter frequency

- Change the following parameter:  
NUC1 [F1] = 31P
- Scroll down to the Program parameter section in the AcqPars window.
- Click the CNST **Edit** button (Constant used in pulse programs).

Program parameters		
L	Edit...	Loop counter
CNST	Edit...	Constant used in pulse programs
CPDPRG	Edit... [Set constants used in pulse programs]	Composite pulse decoupling program (cpd)
PHCOR [degree]	Edit...	Correction angle for phase program
SUBNAM	Edit...	Name of subroutine
ZGOPTNS		Acquisition (zg) options

- Change the following parameter:  
CNST13 = 14 (J 31P/1H long range)

Constant used in pulse programs

CNST[0]	1	CNST[16]	1	CNST[32]	1	CNST[48]	1
CNST[1]	1	CNST[17]	1	CNST[33]	1	CNST[49]	1
CNST[2]	145	CNST[18]	1	CNST[34]	1	CNST[50]	1
CNST[3]	1	CNST[19]	1	CNST[35]	1	CNST[51]	1
CNST[4]	1	CNST[20]	1	CNST[36]	1	CNST[52]	1
CNST[5]	1	CNST[21]	1	CNST[37]	1	CNST[53]	1
CNST[6]	1	CNST[22]	1	CNST[38]	1	CNST[54]	1
CNST[7]	1	CNST[23]	1	CNST[39]	1	CNST[55]	1
CNST[8]	1	CNST[24]	1	CNST[40]	1	CNST[56]	1
CNST[9]	1	CNST[25]	1	CNST[41]	1	CNST[57]	1
CNST[10]	1	CNST[26]	1	CNST[42]	1	CNST[58]	1
CNST[11]	1	CNST[27]	1	CNST[43]	1	CNST[59]	1
CNST[12]	1	CNST[28]	1	CNST[44]	1	CNST[60]	1
CNST[13]	14	CNST[29]	1	CNST[45]	1	CNST[61]	1
CNST[14]	1	CNST[30]	1	CNST[46]	1	CNST[62]	1
CNST[15]	1	CNST[31]	1	CNST[47]	1	CNST[63]	1

Close

The CNST13 long range J value of **14 Hz** is an average value of the two coupling constants **11 Hz** and **17 Hz**, see chapter [Sample ▶ 71](#).

- On the **Spin** button, click the **drop-down** arrow to see more options.
- In the list, select **Turn sample rotation off**.




2D experiments should be run non-spinning.

- On the Workflow button bar, click **Prosol**.

This will load the pulse width and power levels into the parameter set.

- At the command prompt, enter **gradratio**.



## 2D Gradient Experiments


The command executes the AU program **gradratio** to calculate the gradient ratio **GPZ1**, **GPZ2** and **GPZ3** for the nucleus  $^{31}\text{P}$ . To check if the correct gradient ratio values has been entered in the **AcquPars**, follow the steps below.

- In the AcquPars window, scroll down to the Power section.



Parameter	Value	Description
PLW [W]	Edit...	Power level in Watt
PLdB	Edit...	Power level in dB
PLSTRT [dB]	-6	First step for PL switching
PLSTEP	0.1	Step width for PL switching
SHAPE	Edit...	Shaped pulse parameter
GRADIENT	Edit...	Gradient parameters
CAGPARS	Edit...	Parameters for gradient calculation
AMP [%]	Edit...	Amplitude of pulse

- Click **GRADIENT Edit** (Gradient parameters).




Index	GPX [%] (GPX)	GPY [%] (GPY)	GPZ [%] (GPZ)	Filename (GPNAM)
0	0	0	0	
1	0	0	70	SMSQ10.100
2	0	0	30	SMSQ10.100
3	0	0	80.4807	SMSQ10.100

- In the Dataset window, select the **Spectrum** tab.

### 5.8.7 Limit Setting

- On the Workflow button bar, click **SetLimits**.



Close this dialog box after setting frequencies.


1. Open 1D dataset from Browser.
2. Zoom into region of interest.
3. Click OK to set frequencies and return to original dataset.

OK Cancel


- To open the 1D Proton spectrum, right-click on the dataset name in the browser window (e.g. proton exp) and select **Display** or click and hold the left mouse button for dragging the 1D Proton dataset into the spectrum window.

- Expand the spectrum to display all peaks, leaving ca. **1.0 ppm** of baseline on either side of the spectrum.

The solvent peak may be excluded if it falls outside of the region of interest. Digital filtering however is only applied in F2 and the solvent peak will be folding in F1.




- In the setlimits message window, click **OK** to assign the new limit.
- In the message window, click **Close**.




The display changes back to the 2D dataset. Follow the steps below to set the limit in the F1 dimension.

## 2D Gradient Experiments

- On the Workflow button bar, click **SetLimits**.



- To open the 1D  $^{31}\text{P}$  spectrum, right click on the dataset name in the browser window (e.g.  $^{31}\text{P\_exp 1}$ ) and select **Display** or click and hold the left mouse button for dragging the 1D  $^{31}\text{P}$  dataset into the spectrum window.
- Expand the spectrum to display all peaks.



- In the setlimits message window, click **OK** to assign the new limit.
- In the message window, click **Close**.

The display changes back to the 2D dataset.


### 5.8.8 Acquisition

- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.


### 5.8.9 Processing

- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.


This executes a standard processing program **proc2d**. The message shown in the figure above pops up in case of a magnitude 2D experiment and the **apk2d** option is enabled. To disable the **apk2d** option, click the **down arrow** in the **Proc. Spectrum** button in the Workflow button bar and configure the Standard Processing (**proc2d**) program.




- In the apk2 message window, click **Close**.



## 5.8.10 Plotting

Click on any of the 4 icons  to either print the active window, export the active window to a PDF file, copy the active window to the clipboard or, show more publishing options such as E-mailing a dataset, sending a dataset to cloud storage, etc.






## 6 1D Experiments Using Shaped Pulses

### 6.1 Introduction

Selective homonuclear 1D experiments usually start from the selective  $^1\text{H}$  excitation of a given resonance followed by a mixing process. When PFG's are available, the SPFG scheme is highly recommended as a selective excitation scheme. The SPFG or **S**ingle **P**ulsed **F**ield **G**radient **E**cho scheme is a single echo experiment in which the central selective  $180^\circ$  pulse is flanked by two gradient pulses. It is used for efficient selective excitation purposes.




Selective 1D experiments can be easily derived by adding the corresponding mixing process between the SPFG block and the acquisition period.

To run this experiment the instrument has to be equipped with the hardware to do Shaped Pulses and Gradients. Three different ways to run this experiment are discussed in this chapter. Sections 6.3 - 6.5 show how to use the flow bar tools to automatically set up selective excitation regions for selective COSY, NOESY and, TOCSY experiments respectively. Sections 6.6 and 6.7 illustrate how to manually set up selective excitation regions using the on- and off-resonance options.

### 6.2 Sample

A sample of **30mg Menthyl Anthranilate in DMSO-d6** is used for all experiments in this chapter.



## 6.3 1D Selective COSY Using the Flow Bar Tools

### 6.3.1 Introduction


The hard pulses used in all the experiments from the previous chapters are used to uniformly excite the entire spectral width. This chapter introduces soft pulses which selectively excite only one multiplet of a  $^1\text{H}$  spectrum. Important characteristics of a soft pulse include the shape, the amplitude, and the length. The selectivity of a pulse is measured by its ability to excite a certain resonance (or group of resonances) without affecting near neighbors. Since the length of the selective pulse affects its selectivity, the length is selected based on the selectivity desired and then the pulse amplitude (i.e., power level) is adjusted to give a  $90^\circ$  (or  $270^\circ$ ) flip angle.



The transmitter offset frequency of the selective pulse must be set to the frequency of the desired resonance. This transmitter frequency does not have to be the same as  $\omega_{\text{p}}$  (the offset frequency of the hard pulses), but for reasons of simplicity, they are often chosen to be identical.


Most selective excitation experiments rely on phase cycling, and thus subtraction of spectra, to eliminate large unwanted signals. It is important to minimize possible sources of subtraction artifacts, and for this reason it is generally suggested to run selective experiments using pulse field gradients and non-spinning.

Section *1D Selective COSY Experiment Using the On-Resonance Option* [101] describes the acquisition and processing of a one-dimensional  $^1\text{H}$  selective gradient COSY experiment, using the on-resonance option. The standard Bruker parameter set is SELCOGP and includes the pulse sequence **selcogp** shown in the figure below. It consists of the recycling delay, four radio-frequency (RF) pulses and the acquisition time during which the signal is recorded. The first RF pulse is a  $90^\circ$  pulse, followed by a  $180^\circ$  shaped pulse, a  $180^\circ$  hard pulse and finally a  $90^\circ$  pulse. The delay between the  $180^\circ$  and  $90^\circ$  pulse is  $1/4 \cdot J(\text{H,H})$ . The gradient pulses are applied before and after the shaped pulse.



### 6.3.2 Reference Spectrum

Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment, Experiment Setup through Processing*.




### 6.3.3 Selective Excitation Region Set Up



The selective pulse regions are set up using the integration tools. Power and duration of the shape pulses are calculated using the hard 90° pulse in the prosol table.

- On the menu bar, click **Acquire**.
- On the **More** button, click the **drop-down** arrow to see more options.




- In the list, select **Setup Selective 1D Expts.**

The Workflow button bar changes for setting up the 1D selective experiment.




## 1D Experiments Using Shaped Pulses

- On the Workflow button bar, click **1D Selective Experiment Setup**.



This button is only used for the instruction displayed above.


- In the message window, click **Close**.
- Expand the peak at **4.8 ppm**.
- On the Workflow button bar, click **Define Regions**.
- Integrate the multiplet at **4.8 ppm**.







If desired, other peaks can be integrated and a separate dataset will be created for each integral saved in the region file.


- On the Integration toolbar, click **Save/export integration regions** .
- In the list, select **Save the Region to 'reg'**.



- On the toolbar, click **Return do NOT save regions!** .
- In the message window, click **No**.




- On the **Create Dataset** button, click the **drop-down** arrow to see more options.
- In the list, select **Selective gradient COSY**.



## 1D Experiments Using Shaped Pulses


The default parameters are taken from the standard parameter set **SELCOGP**. If desired, the **Gaus1\_180r.1000** pulse can be changed by clicking on the **Change Shape** button in the above window.

- In the SELCOGP window, click **Accept**.



The new dataset is created, and all parameters are automatically set.


- In the sel1d window, click **OK** to start the acquisition.




### 6.3.4 Processing

- On the menu bar, click **Process**.
- On the **Proc Spectrum** button, click the **drop-down** arrow to see more options.
- In the list, select **Configure Standard Processing**.
- Deselect the following options:
  - Auto-Phasing (apk)
  - Set Spectrum Reference (sref)
  - Auto-Baseline correction (abs)
  - Warn if Processed data exist

- In the proc1d window, click **Execute**.



- Manually adjust the phase of the peaks between **3 ppm** and **1 ppm** for an antiphase pattern and if desired the selective peak at **4.8 ppm** can be phased positive.



## 6.3.5 Plotting Two Spectra on the Same Page

- Display the selective COSY spectrum.


- On the toolbar, click **Multiple display**.

The Multiple display toolbar is displayed:




# 1D Experiments Using Shaped Pulses

- Drag the Reference spectrum into the spectral window.



- To adjust the spectra for best fit, use the  $^2s$   $/2s$   $\uparrow s$   $\downarrow s$   $\leftrightarrow s$  toolbar buttons.


Click on any of the 4 icons  to either print the active window, export the active window to a PDF file, copy the active window to the clipboard or, show more publishing options such as E-mailing a dataset, sending a dataset to cloud storage, etc.

## 6.4 1D Selective NOESY Using the Flow Bar Tools

### 6.4.1 Introduction

This experiment consists of three parts:


- **Selective excitation** of the selected resonance using the SPFG block.
- **Mixing period** consisting of the basic  $90^\circ(^1\text{H})$ -delay- $90^\circ(^1\text{H})$  block in phase polarization transfer to other spins via NOE. Purging gradients are usually applied during the mixing period to remove any residual transverse magnetization.
- **Proton detection** as usual.





## 6.4.2 Reference Spectrum


Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment, Experiment Setup through Processing*.



## 6.4.3 Selective Excitation Region Set Up

The selective pulse regions are set up using the integration tools. Power and duration of the shape pulses are calculated using the hard  $90^\circ$  pulse in the prosol table.

- On the menu bar, click **Acquire**.
- On the **More** button, click the **drop-down** arrow to see more options.




- In the list, select **Setup Selective 1D Expts.**

The Workflow button bar changes for setting up the 1D selective experiment.


- On the Workflow button bar, click **1D Selective Experiment Setup**.
- In the message window, click **Close**.

## 1D Experiments Using Shaped Pulses



There is no other function to this button then the instruction displayed above.


- Expand the peak at **4.8 ppm**.
- On the Workflow button bar, click **Define Regions**.





The Define Regions toolbar is displayed:




- Integrate the multiplet at **4.8 ppm**.




If desired, other peaks can be integrated and a separate dataset will be created for each integral saved in the region file.

- On the toolbar, click **Save/export regions** .
- In the list, select **Save Regions to 'reg'**.
- On the toolbar, click **Return do NOT save regions!** .
- In the message window, click **No**.




- On the **Create Dataset** button, click the **drop-down** arrow to see more options.
- In the list, select **Selective gradient NOESY**.



The default parameters are taken from the standard parameter set **SELNOGP**. The mixing time **D8** is dependent on the size of the Molecule and the magnetic strength. It can vary from a large Molecule to a small one from **100 ms** to **800 ms**. If desired, the **Gaus1\_180r.1000** pulse can be changed by clicking on the **Shape** button in the above window.


## 1D Experiments Using Shaped Pulses

- Enter:  
D8 = **0.450**  
NS = **32**
- In the SELNOGP window, click **Accept**.



The new dataset is created, and all parameters are automatically set.


- In the sel1d window, click **OK** to start the acquisition.




### 6.4.4 Processing

- On the menu bar, click **Process**.
- On the **Proc Spectrum** button, click the **drop-down** arrow to see more options.
- In the list, select **Configure Standard Processing**.
- In the proc1d window, deselect the following options:
  - Auto-Phasing (apk)
  - Set Spectrum Reference (sref)
  - Auto-Baseline correction (abs)
  - Warn if Processed data exist

- In the proc1d window, click **Execute**.



- Manually adjust the phase of the selective peak at **4.8 ppm** to show negative absorption to assure the correct phasing of the NOE peaks between **3 ppm** and **1 ppm**. Dependent on the field strength the peaks could be either positive or negative.



## 6.4.5 Plotting Two Spectra on the Same Page


- Display the selective TOCSY spectrum.

- On the toolbar, click **Multiple display**.


The Multiple display toolbar is displayed:



- Drag the Reference spectrum into the spectral window.




- To adjust the spectra for best fit, use the  $\frac{1}{2}$   $\frac{1}{2}$   $\frac{1}{2}$   $\frac{1}{2}$  toolbar buttons.

Click on any of the 4 icons  to either print the active window, export the active window to a PDF file, copy the active window to the clipboard or, show more publishing options such as E-mailing a dataset, sending a dataset to cloud storage, etc.

## 6.5 1D Selective TOCSY Using the Flow Bar Tools


This experiment consists of three parts:

- **Selective excitation** of the selected resonance using the SPFG block.
- **Mixing period** to achieve in phase polarization transfer to other spins. This is usually achieved by applying some isotropic mixing sequence like MLEV, WALTZ or DIPSI pulse trains. This in-phase transfer avoids possible cancellation when the coupling is poorly resolved.
- **Proton detection** as usual.



### 6.5.1 Reference Spectrum


Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment*, *Experiment Setup* through *Processing*.



### 6.5.2 Selective Excitation Region Set Up

The selective pulse regions are set up using the integration tools. Power and duration of the shape pulses are calculated using the hard 90° pulse in the prosol table.

- On the menu bar, click **Acquire**.
- On the **More** button, click the **drop-down** arrow to see more options.




- In the list, select **Setup Selective 1D Expts.**

The Workflow button bar changes for setting up the 1D selective experiment.


## 1D Experiments Using Shaped Pulses

- On the Workflow button bar, click **1D Selective Experiment Setup**.



This button is only used for the instruction displayed above.


- In the message window, click **Close**.
- Expand the peak at **4.8 ppm**.
- On the Workflow button bar, click **Define Regions**.
- Integrate the multiplet at **4.8 ppm**.




If desired, other peaks can be integrated and a separate dataset will be created for each integral saved in the region file.

- On the Integration toolbar, click **Save/export integration regions** .
- In the list, select **Save the Region to 'reg'**.




- On the toolbar, click **Return do NOT save regions!** 
- In the message window, click **No**.



- On the **Create Dataset** button, click the **drop-down** arrow to see more options.
- In the list, select **Selective gradient TOCSY**.

The default parameters are taken from the standard parameter set **SELMLGP**. If desired, the **Gaus1\_180r.1000** pulse can be changed by clicking on the **Shape** button in the above window. A mixing time of **0.06 s** to **0.08 s** is typically for the **TOCSY** experiment.


- Enter:  
D9 = **0.08**  
NS = **8**
- In the SELMLGP window, click **Accept**.



The new dataset is created, and all parameters are automatically set.


## 1D Experiments Using Shaped Pulses

- In the sel1d window, click **OK** to start the acquisition.




### 6.5.3 Processing

- On the menu bar, click **Process**.
- On the **Proc Spectrum** button, click the **drop-down** arrow to see more options.
- In the list, select **Configure Standard Processing**.
- In the proc1d window, deselect the following options:
  - Auto-Phasing (apk)
  - Set Spectrum Reference (sref)
  - Auto-Baseline correction (abs)
  - Warn if Processed data exist
- In the proc1d window, click **Execute**.



- Manually phase all peaks for positive absorption.



## 6.5.4 Plotting Two Spectra on the Same Page


- Display the selective TOCSY spectrum.
- On the toolbar, click **Multiple display**.


The Multiple display toolbar is displayed:




## 1D Experiments Using Shaped Pulses

- Drag the reference spectrum into the spectral window.



- To adjust the spectra for best fit, use the  toolbar buttons.

Click on any of the 4 icons  to either print the active window, export the active window to a PDF file, copy the active window to the clipboard or, show more publishing options such as E-mailing a dataset, sending a dataset to cloud storage, etc.

### 6.5.5 Plotting All 4 Experiments on the Same Page


- Display the selective NOESY spectrum.


- On the toolbar, click **Multiple display**. 

The Multiple display toolbar is displayed:



- Drag the selective COSY spectrum into the spectral window
- Drag the selective TOCSY spectrum into the spectral window
- Drag the Reference spectrum into the spectral window.




Click on any of the 4 icons  to either print the active window, export the active window to a PDF file, copy the active window to the clipboard or, show more publishing options such as E-mailing a dataset, sending a dataset to cloud storage, etc.

## 6.6 1D Selective COSY Experiment Using the On-Resonance Option

### 6.6.1 Reference Spectrum

Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment, Experiment Setup through Processing*.




## 6.6.2 Selective Excitation Region Set Up





Ensure that the SW is large enough to cover the entire spectrum accounting for the position of O1. The shaped pulse is applied on resonance (at the O1 position). The power level and width of the excitation pulse must be known and entered in the Prosol parameter table.


- At the command prompt, type **wrpa**.




- Change NAME = **sel\_cosy**.
- In the wrpa window, click **OK**.
- At the command prompt, type **re** and hit **Enter**.



- Change NAME = **sel\_cosy**.
- In the re window, click **OK**.
- Expand peak at **4.8 ppm**.
- On the toolbar, click **Set RF from cursor**. 




- Move the cursor line into the center of the multiplet.
- To set the frequency, click **left**.
- In the O1/O2/O3 window, click **O1**.




### 6.6.3 Setting Up the Selective COSY

- On the menu bar, click **Start**.
- On the Workflow button bar, click **Read Pars**.
- In the **Find file names** field, enter **SEL\*** to display all selective parameter sets as shown in the figure below.

# 1D Experiments Using Shaped Pulses



- Select **SELCOGP**.
- In the Parameter Sets: rpar window, click **Read**.
- Select the acqu, proc and outd parameter options only.
- In the **Keep parameters** list of values, select **P1, O1, PLW1**.
- Enable the **Keep parameters** option.
- In the rpar window, click **OK**.



- In the Dataset window, select the **Title** tab and enter:  
**1D Selective COSY experiment**  
**30 mg Menthyl Anthranilate in DMSO-d6**
- To store the title, click **Save**.
- In the Dataset window, select the **Spectrum** tab.
- On the menu bar, click **Acquire**.

For the following steps, use the Workflow button bar.

- Click **Spin** and select **Sample rotation off**.





1D selective experiments should be run non-spinning.


- To load the probe/solvent depended parameters, click **Prosol**.

#### 6.6.4 Acquisition


- To start the acquisition, click **Run**.

#### 6.6.5 Processing


- On the menu bar, click **Process**.
- On the **Proc Spectrum** button, click the **drop-down** arrow to see more options.
- In the list, select **Configure Standard Processing**.
- In the proc1d window, deselect the following options:
  - **Auto-Phasing (apk)**
  - **Set Spectrum Reference (sref)**
  - **Auto-Baseline correction (abs)**
  - **Warn if Processed data exist**
- In the proc1d window, click **Execute**.



## 1D Experiments Using Shaped Pulses



- Expand the spectrum from **4 ppm** to **0.5 ppm**.
- To display an antiphase pattern, adjust the **0-order phase** on the peak at **2.0 ppm**.




- To store the phase values, click **Return and Save phased spectrum.** 

Follow the instructions in chapter [Plotting Two Spectra on the Same Page \[ 87 \]](#) to plot two spectra on the same page.

## 6.7 1D Selective NOESY Experiment Using the Off-Resonance Option

### 6.7.1 Reference Spectrum

Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment, Experiment Setup through Processing*.




### 6.7.2 Selective Excitation Region Set Up




This method does not require a large SW. The shaped pulse is applied off resonance (not on the O1 position). The power level and pulse width of the excitation pulse have to be known and entered into the Prosol parameters.

- At the command prompt, type **wrpa**.
- Change NAME = **sel\_noesy**.
- In the wrpa window, click **OK**.




## 1D Experiments Using Shaped Pulses

- At the command prompt, type **re**.



- Change NAME = **sel\_noesy**.
- In the re window, click **OK**.
- Expand the peak at **4.8 ppm**.




- Move the cursor line to the center of the peak.
- **Step 1:** Write down the cursor offset frequency value displayed in the upper left of the spectrum window (e.g. **1439.75**).



To display the cursor information, right-click inside the spectrum window and select **Spectra Display Preferences** and enable **Cursor information** in the Spectra Display Preferences window.

- **Step 2:** At the TopSpin command prompt, type **O1**.



- **Step 3:** Write down the current value (e.g. **1853.43**).
- **Step 4:** Calculate the difference of step 1 and 3 (e.g. **-413.68**).
- In the O1 window, click **Cancel**.



If the signal is down field of O1, a positive value must be entered for spoff. If the signal is up field of O1, spoff will have a negative value.

## 6.7.3 Setting Up the Selective NOESY


- On the menu bar, click **Start**.
- On the Workflow button bar, click **Read Pars**.
- In the **Find file names** field, enter **SEL\*** to display all selective parameter sets as shown in the figure below.



- Select **SELNOGP**.
- In the Parameter Sets: rpar window, click **Read**.
- Select the acqu, proc and outd parameter options only.
- In the **Keep parameters** list of values, select **P1, O1, PLW1**.
- Enable the **Keep parameters** option.

## 1D Experiments Using Shaped Pulses

- In the rpar window, click **OK**.



- In the Dataset window, select the **Title** tab.
- Make the following changes:
  - 1D Selective NOESY experiment**
  - 30 mg Menthyl Anthranilate in DMSO-d6**
- To store the title, click **Save**.
- In the Dataset window, select the **Spectrum** tab.
- On the menu bar, click **Acquire**.

For the following steps, use the Workflow button bar.

- Click **Spin** and select **Sample rotation off**.



1D selective experiments should be run non-spinning.

- To load the probe/solvent depended parameters, click **Prosol**.
- In the Dataset window, select the **AcquPars** tab.
- Make the following changes:
  - PULPROG = **selnogp**
  - D8 = **0.450**
  - DS = **8**
  - NS = **64**
  - SPNAM2 = **Gaus1\_180r.1000**

SPOFF2 = value from *Step 4: Calculate the difference of step 1 and 3* in chapter [Selective Excitation Region Set Up](#) [▶ 109].




The mixing time **D8** is dependent on the size of the molecule and the magnetic strength. It can vary from a large molecule to a small one from **100 ms** to **800 ms**.

## 6.7.4 Acquisition

- To start the acquisition, click **Run**.


## 6.7.5 Processing

- On the menu bar, click **Process**.
- On the **Proc Spectrum** button, click the **drop-down** arrow to see more options.
- In the list, select **Configure Standard Processing**.
- In the proc1d window, deselect the following options:
  - Auto-Phasing (apk)
  - Set Spectrum Reference (sref)
  - Auto-Baseline correction (abs)
  - Warn if Processed data exist
- In the proc1d window, click **Execute**.



- Expand the spectrum from **4 ppm** to **0.5 ppm**.
- Manually adjust the phase of the selective peak at **4.8 ppm** to show negative absorption to assure the correct phasing of the NOE peaks between **3 ppm** and **1 ppm**. Dependent on the field strength the peaks could be either positive or negative.

## 1D Experiments Using Shaped Pulses



- To store the phase values, click **Return and Save phased spectrum.** 

Follow the instructions in chapter [Plotting Two Spectra on the Same Page \[ 87 \]](#) to plot two spectra on the same page.




# 7 2D Experiments using Shaped Pulses

## 7.1 2D Edited HSQC Experiment with Adiabatic Pulses


### 7.1.1 Introduction

The HSQC experiment is the method of choice for a very well resolved H,C correlation. However, in contrast to the HMQC this experiment uses  $^{13}\text{C}$   $180^\circ$  pulses, which causes problems if the  $180^\circ$  pulses are too long in duration (e.g. TXI probes) to cover a very wide spectral range. This leads to phasing problems for high field instruments above 500 MHz. This problem is avoided by applying  $^{13}\text{C}$  frequency-swept adiabatic  $180^\circ$  pulses which can cover the large  $^{13}\text{C}$  spectral width.

The figure below shows the edited HSQC sequence using hard  $^{13}\text{C}$  pulses:




The edited HSQC sequence using shaped pulses for all  $180^\circ$  pulses on the f2-channel with gradients in the back-incept is shown below. For improvement of the phasing the pulse sequence using matched sweep adiabatic pulses **hsqcedetgpsisp2.3** is used in this chapter. This pulse sequence is used in the recommended Bruker parameter set HSQCEDETGPSP-SISP\_ADIA. If desired the sequence **hsqcedetgpsisp2.4** can be used to suppress the COSY peaks.




### 7.1.2 Sample

30 mg Menthyl Anthranilate in DMSO-d<sub>6</sub>



### 7.1.3 Reference spectrum

Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment, Experiment Setup through Processing*.



The reference spectrum is necessary to adjust the spectral limits of the sweep width in the **F2** dimension and to use it for the projection. The HSQCEDETGPSP\_ADIA parameter set has a default sweep width in the **F1** dimension of **165 ppm**. If a Carbon DEPT135 or DEPT45 spectrum of the same sample is available, the F1 sweep width can be further reduced using the **setlimits** AU-program.

### 7.1.4 Setting up the HSQC experiment

The steps below assume that the sample remains in the magnet after observing the proton spectrum.

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the New Dataset window, enter or select:  
NAME = **shape\_hsqc\_exp**  
EXPNO = 1  
Experiment select **HSQCEDETGPSP\_ADIA**  
Set Solvent select **DMSO**
- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.
- To tune the probe, click **Tune**.



The last step is necessary to tune the X-channel which is in this case  $^{13}\text{C}$ . This performs an **atma** (automatic tuning) and requires a probe equipped with an automatic tuning module. Other options can be selected by clicking on the down arrow inside the **Tune** button.


- Click **Spin** and select **Turn sample rotation off**.

2D experiments should be run non-spinning.

- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent depended parameters, click **Prosol**.

### 7.1.4.1 Limit Setting


- On the Workflow button bar, click **SetLimits**.



- To open the 1D Proton spectrum, right-click on the dataset name in the browser window (e.g. proton\_exp 1) and select **Display** or click and hold the left mouse button for dragging the 1D Proton dataset into the spectrum window.
- Expand the spectrum to display all peaks, leaving ca. **1.0 ppm** of baseline on either side of the spectrum.

The solvent peak may be excluded if it falls outside of the region of interest. Digital filtering however is only applied in F2 and the solvent peak will fold in F1.

## 2D Experiments using Shaped Pulses



- In the setlimits message window, click **OK** to assign the new limit.




Figure 7.1:

- In the message window, click **Close**.

The display automatically changes back to the 2D dataset.


### 7.1.4.2 Acquisition

- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.


## 7.1.4.3 Processing

The steps below will guide you through the processing and the manual phase correction on the edited HSQC experiment.

- On the menu bar, click **Process**.
- On the **Proc Spectrum** button, click the **drop-down** arrow to see more options.
- In the list, select **Process F2+F1 Axis (xfb)**, or at the command prompt, type **xfb**.



- On the Workflow button bar, click **Adjust Phase**.


The Adjust phase toolbar is displayed. 

- Select the peak at **7.7 ppm/130.9 ppm**.


## 2D Experiments using Shaped Pulses

- Right-click and on the shortcut menu, select **Add**.


Add  
Remove Row/Col  
Remove All




- Repeat the last step for the peaks at 4.8 ppm/73.2 ppm and 0.76 ppm/16.8 ppm.





- Click **Start the phase correction on rows.** 






- Adjust the **0** order phase on the peak at 7.7 ppm and the **1<sup>st</sup>** order phase on the peak at 0.76 ppm.




## 2D Experiments using Shaped Pulses

- To store the phase values, click **Return and Save phased spectrum.** 
- Click **Start the phase correction on columns.** 
- Adjust the **0** and **1<sup>st</sup>** order phase.



- To store the phase values, click **Return and Save phased spectrum.** 
- To exit the phase window, click the **Return** button. 






## 7.2 2D Selective HMBC Experiment


### 7.2.1 Introduction

The **Semi-selective 2D HMBC experiment** is a simple modification of the 2D HMBC pulse sequence shown in the first figure below in which one of the two carbon  $90^\circ$  pulses is applied selectively on a specified region, see second figure below. The main purpose is to achieve better resolution in the indirect dimension and therefore is recommended when highly overlapped carbon spectra precludes an easy resonance assignment. There are three ways to set this experiment up. Each one will be covered separately below. Before running any one of these methods, you need at least a **Proton** or either a **2D HMBC** or a **1D Carbon spectrum** if possible.



### 7.2.2 Sample


50 mM Cyclosporin in  $C_6D_6$



### 7.2.3 Preparation Experiments

#### 7.2.3.1 1D Proton Experiment


Run a **1D Proton** spectrum of Cyclosporin, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment*, *Experiment Setup* through *Processing* using  $C_6D_6$  as the solvent.



#### 7.2.3.2 2D HMBC Experiment


The steps below assume that the sample remains in the magnet after observing the proton spectrum.

Run a **2D HMBC** experiment of Cyclosporin following the instructions in TopSpin Guide Book *Basic NMR experiments*, chapter *2D HMBC experiment* using  $C_6D_6$  as the solvent.



### 7.2.3.3 1D Proton Decoupled Carbon Experiment

Run a **1D Carbon** spectrum of **Cyclosporin**, following the instructions the TopSpin Guide Book *Basic NMR Experiments*, Chapter *1D Carbon experiment*, Paragraph *Experiment Setup* through *Processing* using  $C_6D_6$  as the solvent.




The carbon spectrum is necessary for method 1 but not for method 2, because the sample concentration is too low to get a  $^{13}C$  spectrum in a reasonable time frame.

### 7.2.3.4 Method 1 for Setting Up the Selective HMBC Experiment

This method requires a 1D Proton decoupled  $^{13}C$  spectrum, if it can be obtained with a reasonable number of scans for adequate S/N (Signal to noise).

- Display the carbon spectrum as observed in the last chapter [1D Proton Decoupled Carbon Experiment](#) [▶ 123].



### 7.2.3.5 Selective Excitation Region Setup

The selective pulse region is set up same way as the 1D selective experiments using the Workflow button bar. Power and duration of the shape pulses are calculated using the hard 90° pulse in the prosol table.

- On the menu bar, click **Acquire**.
- On the **More** button, click the **drop-down** arrow to see more options.



- In the list, select **Setup Selective 1D Expts.**

The Workflow button bar changes for setting up the 1D selective experiment.




- On the Workflow button bar, click **Define Regions**.


The Define Regions toolbar is displayed:



- Integrate the region from **175 ppm** to **167 ppm**.




- On the Integration toolbar, click **Save/export integration regions** .



- In the list, select **Save Regions to reg**.

To exit from the integration mode:


- On the toolbar, click **Return do NOT save regions!** .




- In the message window, click **No**.

## 2D Experiments using Shaped Pulses


- On the **Create Dataset** button, click the **drop-down** arrow to see more options.



- In the list, select **2D Selective HMBC**.
- In the SHMBCCTETGPL2ND window, click **Accept**.



All parameters are automatically calculated and stored as an increment in the next free experiment number of the dataset.



- In the sel1d window, click **OK** to start the acquisition.


**The acquisition starts momentarily.**

## 7.2.3.6 Processing

The pulse program **shmbcctetgpl2nd** is a phase sensitive program. However the data should be processed in magnitude mode. Do not use the **Proc. Spectrum flow** button, rather follow the steps below for the processing.

**Note:** These instructions are at the bottom of the comments in the pulse program file.

- At the command prompt, type **xfb** to process the data in both dimensions.
- At the command prompt, type **xf2m** to calculate magnitude spectrum in F2.





## 7.2.3.7 Method 2 for Setting Up the Selective HMBC Experiment

This method uses a regular **2D HMBC** acquired spectrum for setting up the **2D selective HMBC** experiment. In this example, the **1D Proton decoupled carbon** spectrum is only used to display the F1 projection and is not necessary to obtain the **2D selective HMBC**.

- Display the **HMBC** spectrum as observed in chapter [2D HMBC Experiment](#) [ 122].

## 2D Experiments using Shaped Pulses


- Expand the region including all cross peaks (e.g. 163 ppm to 179 ppm).




- On the menu bar, click **Acquire**.




- On the **More** button, click the **drop-down** arrow to see more options.



- In the list, click the arrow in the selection **Setup Selective 2D Expts.** and on the short-cut menu, select **Band Selective 2D HMBC**.




- In the SHMBCCTETGPL2ND window, click **Accept**.



## 2D Experiments using Shaped Pulses

All parameters are automatically calculated and stored as an increment in the next free experiment number of the dataset.




- To start the acquisition, click **OK**.

### 7.2.3.8 Processing

The pulse program **shmbcctetgpl2nd** is a phase sensitive program however the data should be processed in magnitude mode. Do not use the **Proc. Spectrum flow** button, rather follow the steps below for the processing.

- At the command prompt, type **xfb** to process the data in both dimensions.
- At the command prompt, type **xf2m** to calculate magnitude spectrum in F2.




The **Selective HMBC** has a significantly higher <sup>13</sup>C resolution compared to the standard **HMBC** experiment.

# 8 T1 Experiment


## 8.1 Introduction

The inversion-recovery experiment measures longitudinal or spin-lattice T1 relaxation times of any nucleus.

The basic pulse sequence consists of a 180° pulse that inverts the magnetization to the -z axis. During the following delay, relaxation along the longitudinal plane takes place. Magnetization comes back to the original equilibrium z-magnetization. A 90° pulse creates transverse magnetization. The experiment is repeated for a series of delay values taken from a variable delay list. A 1D spectrum is obtained for each value of vd and stored in a pseudo 2D dataset. The relaxation time d1 must be set to 5\*T1. A rough estimation of the T1 value can be calculated from the null-point value by using  $T1 = t_{null} / \ln(2)$ .



## 8.2 Sample




Menthyl Anthranilate in DMSO-d<sub>6</sub>

## 8.3 Proton Inversion-Recovery T1 Experiment

### 8.3.1 Preparation Experiment

Run a 1D Proton spectrum, following the instructions in the *TopSpin Guide Book Basic NMR Experiments*, Chapter *1D Proton experiment*, Paragraph *Experiment Setup* through *Processing*.



The reference spectrum is necessary to adjust the spectral limits of the sweep width to gain more data points.

## 8.3.2 Setting up the T1 Experiment

The steps below assume that the sample remains in the magnet after observing the proton spectrum.

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the New Dataset window, enter or select:  
NAME = **t1\_exp**  
EXPNO = **1**  
Experiment: select **PROTONT1**  
Set Solvent: select **DMSO**
- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.

For the following steps, use the Workflow button bar.


- To tune the probe, click **Tune**.
- Click **Spin** and select **Sample rotation off**.

T1 experiments should be run non-spinning.

- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent depended parameters, click **Prosol**.

### 8.3.3 Limit Setting

- On the Workflow button bar, click **SetLimits**.




To open the 1D Proton spectrum

- Right click on the dataset name in the browser window (e.g. **proton\_exp 1**) and select **Display**


or

- Click and hold the left mouse button and drag the 1D Proton dataset into the spectrum window.
- Expand the spectrum to display all peaks, leaving about **1.0 ppm** of baseline on either side of the spectrum.

The solvent peak may be excluded if it falls outside of the region of interest.




- Click **OK** in the setlimits message window to assign the new limit.




- In the message window, click **Close**.

The display changes back to the 2D dataset.

- In the Dataset window, select the **AcquPars** tab.
- Click **Show pulse program parameters**. 
- Make the following changes:
  - D1 = **15**
  - VDLIST = **t1delay\_MA**



- Click **Edit variable delay list** right of the VDLIST name box.




- Enter the variable delay values as shown in the figure above.
- Click **File** and **Save**.
- Click **File** and **Close**.
- In the Dataset window, select the **Spectrum** tab.


### 8.3.4 Acquisition

- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

### 8.3.5 Processing


- On the menu bar, click **Process**.
- At the command prompt, type **rser 10**.
- At the command prompt, type **ef**.
- On the Workflow button bar, click **Adjust Phase**.
- Adjust the phase manually or enter **apk** for automatic phase correction.




- On the Adjust Phase toolbar, click **Save for nD spectrum**. 

- On the toolbar, click **Return, do NOT save phased spectrum**. 


The spectrum will go back to the un-phased view since the phase correction values were stored only for the 2D spectrum.

- On the toolbar, click **To Last 2D data** to go back to the 2-D spectrum display. 
- At the command prompt, type **xf2** to process only the F2 axis.
- Type **abs2** to baseline correct the rows.




## 8.3.6 T1 Calculation

- On the menu bar, click **Applications**.
- On the **Dynamics** button, click the drop-down arrow to see more options and in the list, select **T1/T2 Module**.




The flow buttons change to determine the T1 / T2 relaxation times:




While executing the steps below, message windows will be displayed. Please read each message thoroughly and follow the instructions.

- On the Workflow button bar, click **Fid**.
- In the Extract a row from 2d data window, click **Spectrum**.







- Enter Slice Number = **10**.




- In the message window, click **OK**.




 **Peaks/Ranges**

- On the Workflow button bar, click **Peaks/Ranges**.
- In the Define Peaks and/or Integrals window, click **Manual Integration**.




# T1 Experiment


- In the Prepare relaxation data window, click **OK**.




- Define the regions by drawing an integral over the peaks of interest.



- On the Integration toolbar, click **Save/export integration regions** .
- In the list, select **Export Region To Relaxation Module**.




- On the Workflow button bar, select **Relaxation**.



By default, the selected areas are peak-picked, and the first peak is displayed in the Relaxation window.

- On the Workflow button bar, select **Fitting**.



- In the message window, click **Close**.

**Relaxation parameters**

**General Parameters**

16	FID # for phase determination
1000.0	Left limit for baseline correction
-1000.0	Right limit for baseline correction
5	Number of drift points
1.0E-5	Convergence limit
16	Number of points
1	First slice
1	Slice increment
1.0	Peak sensitivity

**Fitting Function**

uxnmrt1	Function Type
1	Number of components
vclist	List file name
0.001	Increment (auto)
pd	to pick data points

**Iteration control parameters**


Guesses    Reset

**Additional Parameters**



10000.0	GAMMA(Hz/G)
10.0	LITDEL(msec)
100.0	BIGDEL(msec)
1.0	GRADIEN(G/cm)

OK    Apply    Cancel

- In the Relaxation parameters window, click **OK** and select **Area** as Fitting type.



- On the Workflow button bar, select **Calculation**. > Calculation
- In the message window, click **Close**.



- In the T1/T2 tools bar, click **Calculate fit for all peaks.** 

**Brief Report**

Region 1 from 7.797 to 7.650 ppm  
T1 = 1.339s

Region 2 from 7.313 to 7.159 ppm  
T1 = 1.294s


Region 3 from 6.860 to 6.413 ppm  
T1 = 555.498m


Region 4 from 4.891 to 4.725 ppm  
T1 = 644.916m

Region 5 from 3.414 to 3.321 ppm  
T1 = 1.110s

Region 6 from 2.076 to 1.372 ppm  
T1 = 398.815m

Region 7 from 1.171 to 0.552 ppm  
T1 = 378.896m

- On the Workflow button bar, select **Report.** 



```

1 Dataset :
2 C:/Bruker/TopSpin/examdata/t1_exp/1/pdata/1
3 AREA fit :
4 I[t]=I[0]+P*exp(-t/T1)
5
6 16 points for Integral 1, Integral Region from 7.797 to 7.650 ppm
7 Results Comp. 1
8
9 I[0] = 1.047e+00
10 p = -2.003e+00
11 T1 = 1.339s
12 SD = 7.030e-04
13
14 tau ppm integral intensity
15
16 1.000m 7.698 -5.4489e+08 -6.8373e+07
17 50.000m 7.698 -5.0393e+08 -6.3011e+07
18 100.000m 7.698 -4.6377e+08 -5.7767e+07
19 200.000m 7.698 -3.8748e+08 -4.7869e+07
20 300.000m 7.698 -3.1673e+08 -3.8814e+07
21 500.000m 7.698 -1.8933e+08 -2.298e+07
22 800.000m 7.698 -3.0687e+07 -3.7791e+06
23 1.000s 7.698 5.6667e+07 6.6579e+06
24 1.500s 7.698 2.2514e+08 2.658e+07
25 2.000s 7.698 3.4075e+08 4.0165e+07
26 2.500s 7.698 4.2051e+08 4.9459e+07
27 3.000s 7.698 4.7569e+08 5.5909e+07
28 3.500s 7.698 5.1425e+08 6.042e+07
29 4.000s 7.698 5.406e+08 6.3606e+07
30 4.500s 7.698 5.5853e+08 6.5961e+07
31 5.000s 7.698 5.7096e+08 6.7579e+07
32
  
```



## 9 Pulse Calibration

### 9.1 Introduction



This chapter describes the pulse calibration procedures for determining the  $90^\circ$  transmitter pulse of  $^1\text{H}$ ,  $^{13}\text{C}$  and  $^{15}\text{N}$  nuclei.




Note: If your system has been cortabed, it is always a good practice to obtain spectra with the power check turned on.

### 9.2 Sample

Mixture **0.1 M** each of  $^{15}\text{N}$  enriched Urea and  $^{13}\text{C}$  enriched methanol in **DMSO- $d_6$**  (NMR pulse calibration reference standard at <https://bruker-labscape.store>)




### 9.3 $^1\text{H}$ $90^\circ$ Transmitter Pulse



The pulse program **zg** is used to determine the  $^1\text{H}$   $90^\circ$  transmitter pulse. The sequence consists of one channel **f1** with a recycle delay **d1**, a  $^1\text{H}$  pulse **p1**, followed by the  $^1\text{H}$  signal detection. The signal has maximum intensity if **p1** is a  $90^\circ$  pulse and 2 nulls at a  $180^\circ$  and  $360^\circ$  pulse. A methanol signal region from **3.5 ppm** to **2.8 ppm** is used for this experiment.


#### 9.3.1 Preparation Experiment

Run a **1D Proton** spectrum of urea/methanol in DMSO- $d_6$ , following the instructions from the *TopSpin Guide Book Basic NMR Experiments*, chapter *1D Proton Experiment*, paragraphs *Experiment Setup* through *Processing*.




## 9.3.2 Parameter Setup

- At the command prompt, type **wrpa**.



- In the field *Name*, enter **p90\_proton**.
- Click **OK**.
- At the command prompt, type **re**.







- In the field *Name*, enter **p90\_proton**.
- Click **OK**.

Normally a single on-resonance peak is used to determine the 90° transmitter pulse. For practical reasons the methanol signal region from **3.5 ppm** to **2.8 ppm** is used to measure the <sup>1</sup>H 90° transmitter pulse, since the same signals will also be used in determining the <sup>13</sup>C 90° decoupler pulse.

- Expand the spectrum for the region between **3.5 ppm** and **2.8 ppm**.



- On the toolbar, click **Set sw to current region** and **SFO1 to center of region**. 



- In the pop up window, click **Close**.
- In the Dataset window, select the **AcquPars** tab.
- Enter:  
PULPROG = **zg**

TD = 4096

NS = 1

DS = 0

D1 = 10

- In the Dataset window, select the **ProcPars** tab.
- Enter:  
SI = 2048  
PH\_mod = pk
- In the Dataset window, select the **Spectrum** tab.

## 9.3.3 Acquisition

- On the menu bar, click **Acquire**.
- On the Workflow button bar, click **Gain**.


Alternatively type **rga** at the TopSpin command prompt. To adjust the receiver gain manually, click the **drop-down** arrow on the **Gain** button.

- On the Workflow button bar, click **Run**.

Alternatively, type **go** at the TopSpin command prompt. On the **Go** button, click the **drop-down** arrow to see more options.

## 9.3.4 Processing


- Process and phase correct the spectrum.
- Display the full spectrum.



- Right-click in the spectrum window.
- In the list, select **Save Display Region To...**



- Select **Parameters F1/2 (e.g. used by restore display, ...) [dpl]**.




- Click **OK**.
- At the command prompt, type **wpar H1p90\_urea all** to store the parameter set for future use.

### 9.3.5 Determine the $^1\text{H}$ 90° Transmitter Pulse


- At the command prompt, type **popt** to display the Parameter **OPT**imization window.
- Enter or select from the list boxes:
  - OPTIMIZE = **Step by step**
  - PARAMETER = **p1**
  - OPTIMUM = **POSMAX**
  - STARTVA = **2**
  - NEXP = **20**
  - VARMOD = **LIN**
  - INC = **2**

- Click **Save**.




The **ENDVAL** parameter has been updated.

- In the popt window, click **Start optimize** to display the poptau window.
- Enter **y**.




- Click **OK**.

The parameter optimization starts. The spectrometer acquires and processes 20 spectra by incrementing the parameter p1 from **2  $\mu$ s** by **2  $\mu$ s** to a final value of **40  $\mu$ s**. For each of the 20 spectra, only the spectral region defined above is plotted, and all the spectra are plotted side-by-side in the file *p90\_proton/1/pdata/999* as shown:



- In the Dataset window, select the **Title** tab.




The POSMAX value of p1 is displayed in the title tab window which is the 90° pulse, along with the experiment number and the NEXP value. Write this value down. To obtain a more accurate 90° pulse measurement, follow the steps below:

- Close the popt setup window.
- At the command prompt, type **re 1 1**
- At the command prompt, type **p1**
- Enter the value which corresponds to a 360° pulse (the second zero crossing in the popt spectrum, which should be approximately 4 times the POSTMAX value).
- Step 1: At the command prompt, type **zg** to start the acquisition.
- Step 2: At the command prompt, type **efp**
- Change p1 slightly and repeat steps 1 and 2, until the signal undergoes a zero crossing as expected for an exact 360° pulse.

The signals are negative for a pulse angle slightly less than 360° and positive when the pulse angle is slightly more than 360°.

- Divide the determined 360° pulse value by 4. This will be the exact 90° pulse length for the proton transmitter on the current probe.


## 9.4 $^{13}\text{C}$ 90° Decoupler Pulse



The pulse program used in this procedure is the **decp90** sequence shown in the figure above. The sequence consists of two channels f1 (I) and f2 (S), where in this case f1 is set for  $^1\text{H}$  and f2 to  $^{13}\text{C}$ . Channel f1 shows a recycle delay **d1** followed by a  $90^\circ$  pulse and a delay **d2** =  $1/(2J_{\text{XH}})$  for the creation of antiphase magnetization. A  $^{13}\text{C}$  pulse on channel f2 is been executed after the delay **d2** and then the  $^1\text{H}$  signal is detected. When the  $^{13}\text{C}$  pulse is exactly  $90^\circ$ , the  $^1\text{H}$  signals will go through a null. The methanol signal region from **3.5 ppm to 2.8 ppm** is used for this experiment.

### 9.4.1 Preparation Experiment


Run a **1D** proton spectrum of urea/methanol in DMSO-d<sub>6</sub>, following the instructions the *TopSpin Guide Book Basic NMR Experiments, Chapter 1D Proton Experiment, Experiment Setup* through *Processing*.




### 9.4.2 Parameter Setup

- At the command prompt, type **wrpa** and press **Enter**.


- Change NAME = **p90\_carbon**.
- In the wrpa window, click **OK**.




- At the command prompt, type **re** and press **Enter**.
- Change NAME = **p90\_carbon**
- In the re window, click **OK**.



- Expand the spectrum for the region between **3.5 ppm** and **2.8 ppm**.



- On the toolbar, click **Set sw to current region** and **SFO1 to center of region**.
- In the message window, click **Close**.




- In the Dataset window, select the **AcquPars** tab.
- Enter:  
 PULPROG = **dec90**  
 TD = **2048**  
 NS = **1**  
 DS = **0**
- In the Nucleus2 section of the AcquPars, click **Edit** next to NUC2.

Nucleus 2			
NUC2	off	Edit...	2nd nucleus
O2 [Hz]	1853.43		Frequency offset of 2nd nucleus
O2P [ppm]	6.175		Frequency offset of 2nd nucleus
SFO2 [MHz]	300.1318534		Frequency of 2nd nucleus
BF2 [MHz]	300.1300000		Basic frequency of 2nd nucleus

- Select **<sup>13</sup>C** for **NUC2**.



- In the Edit Spectrometer Parameter window, click **Default** to set the routing.




- In the Edit Spectrometer Parameter window, click **Save**.
- In the **AcquPars** enter:
  - O2[ppm] = **49.5**
  - D1 = **10**
  - CNST2 = **139**
  - P3 = **3**
- On the menu bar, click **Acquire**.
- On the Workflow button bar, click **Prosol**.
- In the Dataset window, select the **ProcPars** tab.
- Enter:
  - SI = **2048**
- In the Dataset window, select the **Spectrum** tab.
- At the command prompt, type **wpar C13p90\_urea all** to store the parameter set for future use.

### 9.4.3 Determine the $^{13}\text{C}$ 90° Decoupler Pulse


- On the menu bar, click **Acquire**.
- On the Workflow button bar, click **Tune**.
- On the Workflow button bar, click **Run**.
- Process and phase correct the spectrum.

# Pulse Calibration




Phase the left doublet negative and the right doublet positive. The water peak at **3.3 ppm** can be ignored and does not have to be in phase.

- Increase **p3** in increments of **1** or **2 μs**, execute **zg** followed by the command **efp** until the signals go through a null or a phase change. This will be the  $^{13}\text{C}$  90° decoupler pulse.




## 9.5 $^{15}\text{N}$ 90° Decoupler Pulse




The pulse program used in this procedure is the **decp90** sequence shown in the figure above. The sequence consists of two channels f1 (I) and f2 (S), where in this case f1 is set for  $^1\text{H}$  and f2 to  $^{15}\text{N}$ . Channel f1 shows a recycle delay  $d1$  followed by a  $90^\circ$  pulse and a delay  $d2 = 1/(2J_{\text{XH}})$  for the creation of antiphase absorption. A  $^{15}\text{N}$  pulse on channel f2 is executed after the delay  $d2$  and then the  $^1\text{H}$  signal is detected. When the  $^{15}\text{N}$  pulse is exactly  $90^\circ$ , the  $^1\text{H}$  signals will go through a null. The urea signal region from **5.6 ppm** to **5.1 ppm** is used for this experiment.

If your system is equipped with a 3rd channel for  $^{15}\text{N}$  observation, you can still follow the same instructions in this chapter with the exceptions of using the pulse sequence **decp90f3** shown in the figure below and the routing which is illustrated in the next section.




### 9.5.1 Parameter Setup

- At the command prompt, type **wrpa** and press **Enter**.
- Change NAME = **p90\_nitrogen**
- In the wrpa window, click **OK**.




## Pulse Calibration


- At the command prompt, type **re** and press **Enter**.
- Change NAME = **p90\_nitrogen**.
- In the re window, click **OK**.



- Expand the spectrum for the region between **5.6 ppm** and **5.1 ppm**.



- On the toolbar, click **Set sw to current region** and **SFO1 to center of region**. 




- In the message window, click **Close**.
- In the Dataset window, select the **AcquPars** tab.
- Enter:  
 PULPROG = **decpg90**  
 TD = **2048**  
 NS = **1**  
 DS = **0**
- In the **Nucleus2** section of the AcquPars, click **Edit** next to **NUC2**.


▼ Nucleus 2		
NUC2	off	Edit...
O2 [Hz]	1853.43	2nd nucleus
O2P [ppm]	6.175	Frequency offset of 2nd nucleus
SFO2 [MHz]	300.1318534	Frequency offset of 2nd nucleus
BF2 [MHz]	300.1300000	Frequency of 2nd nucleus
		Basic frequency of 2nd nucleus

## 9.5.1.1 Two Channel System

- Select **<sup>15</sup>N** for **NUC2**.




- In the Edit Spectrometer Parameter window, click **Default** to set the routing.




- In the Edit Spectrometer Parameter window, click **Save**.

## 9.5.1.2 Three Channel System

- Select <sup>15</sup>N for **NUC2**.



- In the Edit Spectrometer Parameter window, click **Default** to set the routing.




- In the Edit Spectrometer Parameter window, click **Save**.
- In the AcqPars make the following change:  
 $O2[\text{ppm}] = 76$   
 $D1 = 10$   
 $CNST2 = 88.5$   
 $P3 = 6$
- On the Workflow button bar, click **Prosol**.
- In the Dataset window, select the **ProcPars** tab.
- Make the following change:  
 $SI = 2048$
- In the Dataset window, select the **Spectrum** tab.
- At the command prompt, type `wpar N15p90_urea all` to store the parameter set for future use.

## 9.5.2 Determine the $^{15}\text{N}$ 90° Decoupler Pulse


- On the Workflow button bar, click **Tune**.
- On the Workflow button bar, click **Run**.
- Process and phase correct the spectrum.

# Pulse Calibration

- Phase the left side signal negative and the right-side signal positive!



- Increase **p3** in increments of 1 or 2  $\mu\text{s}$ , execute **zg** followed by the command **efp** until the signals go through a null or a phase change. This will be the  $^{15}\text{N}$   $90^\circ$  decoupler pulse.






# 10 <sup>1</sup>H Homonuclear Decoupling

## 10.1 Sample

---

30mg Menthyl Anthranilate in DMSO-d6



## 10.2 <sup>1</sup>H Homonuclear Decoupling Experiment

---


### 10.2.1 Introduction

---

The homonuclear decoupling (homodecoupling) simplifies multiplet structures by irradiating a specific <sup>1</sup>H resonance. Unambiguous assignments and measurement of <sup>1</sup>H-<sup>1</sup>H coupling constants can be performed by analyzing the resulting residual multiplets.

During an homo-decoupling experiment, a conventional <sup>1</sup>H spectrum is recorded. From a second channel, low-power irradiation is applied on a predefined frequency during the acquisition period.

Important parameters to consider are the offset and the power level of the irradiation. It is useful to have a calibration of the field strength delivered from the decoupler in order to optimize the required selectivity and to minimize Bloch-Siegert shift effects.



### 10.2.2 Acquisition

---




The selective pulse regions are set up using the integration tools. Power and duration of the shape pulses are calculated using the hard 90° pulse in the prosol table.

---

# <sup>1</sup>H Homonuclear Decoupling

The steps below assume that the sample remains in the magnet after acquiring the proton experiment.

- On the menu bar, click **Acquire**.
- On the **More** button, click the **drop-down** arrow to see more options.




- In the list, select **Setup Selective 1D Expts.**

The Workflow button bar changes for setting up the 1D selective experiment.




- On the Workflow button bar, click **1D Selective Experiment Setup**.



This button is only used for the instruction displayed above.


- In the message window, click **Close**.
- Expand the peak at **7.7 ppm**.
- On the Workflow button bar, click **Define Regions**.


- Integrate the multiplet at 7.7 ppm.




If desired, other peaks can be integrated and a separate dataset will be created for each integral saved in the region file.

- On the Integration toolbar, click **Save/export integration regions** .
- In the list, select **Save the Region to 'reg'**.




- On the toolbar, click **Return do NOT save regions!** .
- In the message window, click **No**.




# **<sup>1</sup>H Homonuclear Decoupling**

- On the **Create Dataset** button, click the **drop-down** arrow to see more options.
- In the list, select **1H Homonuclear Decoupling**.




- In the PROHOMODEC window, click **Accept**.



The new dataset is created, and all parameters are automatically set.

- In the sel1d window, click **OK** to start the acquisition.




## **10.2.3 Processing**

- On the menu bar, click **Process**.
- On the Workflow button bar, click **Proc Spectrum**.

This executes a processing program including commands such as an exponential window function em, Fourier transformation ft, an automatic phase correction apk and a baseline correction abs.

To configure the commands, click the **drop-down** arrow on the **Proc Spectrum** button and select **Configure Standard Processing**.

- Expand the region from 8 ppm to 6.2 ppm.



The multiplet at **6.4 ppm** should collapse from a triplet to a doublet. If the triplet is partially collapsed, increase the decoupling power pl24 and repeat the steps in chapter Acquisition and [Processing \[ 164\]](#).

## NOTICE

### Material Damage Due to Excessive Power

The NMR probe can be severely damaged if too much power or power over a too long time is applied.

- ▶ Always start to optimize pulses with low power values and short pulses. Respect the pulse and power limits as programmed into the PICS data of the probe.

## 10.2.4 Plotting Two Spectra on the Same Page

- Display the decoupled spectrum.


# <sup>1</sup>H Homonuclear Decoupling

- On the toolbar, click **Multiple display**. 


The Multiple display toolbar is displayed:



- Drag the Reference spectrum (1D proton) into the spectral window.



- To adjust the spectra for best fit, use the  toolbar buttons.

Click on any of the 4 icons  to either print the active window, export the active window to a PDF file, copy the active window to the clipboard or, show more publishing options such as E-mailing a dataset, sending a dataset to cloud storage, etc.

# 11 Proton DOSY Experiment

## 11.1 Introduction

The **DOSY (Diffusion-Ordered Spectroscopy) experiment** provides accurate, noninvasive, molecular diffusion measurements on biofluids, complex chemical mixtures and multi component solutions. In DOSY spectra, chemical shift is along the detected F2 axis and diffusion coefficient is along the other F1 axis.

Molecules in the solution state move. This translational motion is known as Brownian molecular motion and is often simply called diffusion or self-diffusion. Molecular diffusion depends on a lot of physical parameters like size and shape of the molecule, temperature and viscosity.

Pulsed field gradient NMR spectroscopy can be used to measure translational diffusion. By use of a gradient pulse, molecules can be spatially labeled. After this encoding gradient pulse ( $\delta$ ), molecules move during the diffusion time ( $\Delta$ ). Their new position can be decoded by a second gradient pulse. This encoding/decoding procedure results in an attenuation of the NMR signal which can be described by the following equation:

$$I(g) = I(o) \exp \left[ -(\gamma g \delta)^2 D \left( \Delta - \frac{\delta}{3} \right) \right]$$

Where **I** is the observed intensity, **D** is the diffusion coefficient,  $\gamma$  is the gyro magnetic ratio of the encoded nucleus, **g** is the gradient strength,  $\delta$  is the length of the gradient pulse, and  $\Delta$  as mentioned previously is the diffusion time.

The diffusion experiment described below records a series of 1D  $^1\text{H}$  spectra at increasing gradient strengths ( $g$ ) and then fits the signal intensity decay to the above equation to obtain **D**.

## 11.2 Sample

Mixture of Ibuprofen, Pamoic acid and Pinene in DMSO-d6.




The experimental parameters of  $\delta$  (p130) and  $\Delta$  (d20) described here are for this sample. If using a different sample, they will likely be different.

## 11.3 DOSY Experiment


### 11.3.1 Pulse Programs

The DOSY pulse program used in this chapter is a Stimulated spin-echo experiment using bipolar gradients and an additional delay just prior to detection for the ring-down of any possible eddy currents (led).

The figure below is a 1D version of the pulse program and is used to optimize parameters, see the chapter [Parameter Setup \[ 171 \]](#) and.



The pulse program in the figure below is used for the DOSY experiment, see the chapter [Running the Experiment \[ 175\]](#). The difference between these 2 pulse sequences is that the one shown in the figure below is a pseudo 2D sequence and includes the code to automatically increment the gradient strength.



To run this experiment the instrument has to be equipped with the hardware to run gradient experiments.

## 11.3.2 Preparation Experiment

- On the menu bar, click **Acquire | Create Dataset** to open the Create New Dataset window.
- In the New Dataset window, enter or select:
  - NAME = **DOSY\_exp**
  - EXPNO = **1**
  - Experiment: select **PROTON**
  - Set Solvent: select **DMSO**



## DIR

The directory (DIR) is specific to how the data are stored and therefore may show different entries as the one in the figure above. Click the drop-down arrow to browse for a specific directory.

## Title

In the TITLE window enter a text stating the experiment, sample, the solvent and any other useful information. The title information can be used to search for a dataset.

- In the New Dataset window, click **OK**.
- On the menu bar, click **Acquire**.

**For the following steps, use the Workflow button bar.**

- Click **Sample** and eject the sample, if there is one inserted, and insert the new sample.
- Click **Lock** and select **DMSO** solvent.
- To tune the probe, click **Tune**.
- Click **Spin** and select **Turn sample rotation off**.



DOSY experiments should be run non-spinning.

- To autoshim the sample with TopShim for best homogeneity, click **Shim**.
- To load the probe/solvent depended parameters, click **Prosol**.

## 11.3.3 Acquisition

---

- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

## 11.3.4 Processing

---

- Process and phase correct the spectrum.


## 11.3.5 Limit Settings


---




Changing the sweep width to a smaller value increases the resolution.

---


- On the command line, type:  
**wrpa 2**  
**re 2**
- Expand the spectrum from **9 ppm** to **0 ppm**.
- On the toolbar, click **Set sw to current region** and **SFO1 to center of region**. 
- Click **Close**.





### 11.3.6 Parameter Setup


For an accurate DOSY experiment, certain parameters need to be calibrated for each sample to ensure that the observed signal decay is appropriate. This section will walk you through this process.

- On the command line, type **iexpno**.
- In the Dataset window, select the **AcquPars** tab.
- Click **Show pulse program parameters**. 
- Enter:
 

```
PULPROG = ledbpgp2s1d
D20[s] = 0.1
D21[s] = 0.005
GPNAM6 = SMSQ10.100
GPNAM7 = SMSQ10.100
GPNAM8 = SMSQ10.100
GPZ6[%] = 2
GPZ7[%] = -17.13
GPZ8[%] = -13.17
P30[us] = 1400
```
- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.
- Process and phase correct the spectrum.


# Proton DOSY Experiment

- Expand the spectrum to display all peaks.




The asymmetry of the peaks is an artifact which is caused by non-optimal lock settings. The problem can be fixed by applying a loop adjust of the lock especially that of the lock phase.

- At the command prompt, type **loopadj** and confirm with **OK**.
- To start the acquisition, click **Run**.
- Process and phase correct the spectrum.



- Display the full spectrum.

This is the first spectrum (2%). Because the DOSY experiment is fitting the signal as it decays, it is crucial that there be enough signal in this first experiment that it can decay by a factor of about 95%, but still have good enough signal to noise so that these attenuated values are still reasonably error free. Thus, it is recommended that this first experiment have a s/n of approximately 100:1 or greater. If this is not the case, increase the NS.




The steps described below are necessary to make sure the attenuation level is sufficient at the final (95%) gradient strength. In this second (95% gradient strength) experiment, there should still be signal, but it should be attenuated by a factor of 90-95% as compared to that of the first (2% gradient strength). If there is only noise, then p30 and or d20 need to be reduced. If there is less than 90-95% attenuation, p30 and or d20 need to be increased.

For both the cases of increasing or decreasing the gradient pulse/delays adjusting p30 will have more of an effect than adjusting d20.


Keep in mind that the recommended safety limit for p30 is 2.5 ms, after this limit has been reached, further attenuation must be achieved through increasing d20.

If large changes to p30 and d20 are necessary, then it is recommended to re-run both experiments because these values will affect the intensity of not only the attenuation in the 2nd spectrum, but the starting intensity in the 1st spectrum.

- On the command line, type **iexpno**.
- In the Dataset window, select the **AcquPars** tab.
- Click **Show pulse program parameters**. 

# Proton DOSY Experiment

- Enter:  
GPZ6[%] = **95**
- To start the acquisition, click **Run**.
- Process and phase correct the spectrum.




If there are no signals present, then p30 and/or d20 need to be reduced.

- On the toolbar, click **Multiple display**. 

The Multiple display toolbar is displayed:



- Drag the previous experiment into the multiple display window (in this example it is experiment # 3) or type **re 3**.




As described above, there need to be signal remaining in the 2nd experiment, but the intensity difference of the two spectra should be a factor of 90-95%. If neither of these are true, it is necessary to change p30, and d20 accordingly.

- To exit the multiple display, click **Return**. ↵


## 11.3.7 Running the Experiment

- On the command line, type **ixpno**.
- In the Dataset window, select the **AcquPars** tab.
- Enter:  
PULPROG = **ledbpgp2s**
- Click **Change acquisition dimension of current dataset**. ↓ 1,2.




## Proton DOSY Experiment

- In the parmmod window, select **2D** and click **OK**.
- Enter:  
TD[F1] = **25**  
FnMODE = **QF**
- At the command prompt, type **dosy**.



- In the dosy window, click **Close**.




- Enter **2** for first gradient amplitude and click **OK**.




- In the dosy window, enter **95** for final gradient amplitude and click **OK**.






- In the dosy window, enter **25** for the number of points and click **OK**.




- In the dosy window, enter **l** for the ramp type and click **OK**.



- To start the acquisition, click **OK**.



## 11.3.8 Processing



- In the Dataset window, select the **Fid** tab.



This step is only used to illustrate the DOSY experiment as a decay function.

- In the Dataset window, select the **ProcPars** tab.
- Make the following changes:
  - SI [F1] = **256**
  - PH\_mod [F1] = **no**
  - PH\_mod [F2] = **pk**
- At the command prompt, type **rser 1** to read in the first serial file of the 2D experiment.
- At the command prompt, type **em** to apply the window function.
- At the command prompt, type **ft**.
- On the menu bar, click **Process**.
- On the Workflow button bar, click **Adjust Phase**.
- Process and phase correct the spectrum.
- On the Adjust Phase toolbar, click **Save for nD spectrum.** 
- On the toolbar, click **Return, do NOT save phased spectrum.** 




The spectrum will go back to the un-phased view since the phase correction values were stored only for the 2D spectrum.

- On the toolbar, click **To Last 2D data** to go back to the 2-D spectrum display. **2D**
- At the command prompt, type **xf2**.
- At the command prompt, type **abs2**.
- At the command prompt, type **setdiffparm**.



This command transfers experimental parameters into the values used for fitting the data.




### 11.3.9 Calculating the Diffusion Coefficient


As you follow the steps below, message windows with important instructions will pop up. Please read these instructions very carefully.

- On the menu bar, click **Applications**.

- On the **Dynamics** button, click the **drop-down** arrow to see more options.
- In the list, select **T1/T2**.



The Workflow buttons change to the mode: Determination of the **T1 / T2** relaxation times.




- On the Workflow button bar, click **Fid**.





While executing the next steps, message windows will pop up. Please read each message thoroughly and follow the instructions in it.

- In the Extract a row from 2d data window, click **Spectrum**.




- In the field Slice Number, enter **1** and in the message window click **OK**.






- On the Workflow button bar, click **Peaks/Ranges**.
- In the Define Peaks and/or Integrals window, click **Manual Integration**.





- In the Prepare relaxation data window, click **OK**.




# Proton DOSY Experiment


- Define the regions by drawing an integral over the peaks of interest.




- On the Integration toolbar, click **Save/export integration regions** .
- From the drop-down list, select **Export Region To Relaxation Module**.
- On the Workflow button bar, click **Relaxation**.
- In the group box Fitting type, select **Area**.




- On the Workflow button bar, click **Fitting**.
- In the message window, click **Close**.




- In the Fitting Function group box, select **vargrad** and **difflist** and click **OK**.



- On the Workflow button bar, click **Calculation**.
- In the message window, click **Close**.



- On the T1/T2 toolbar, click **Calculate fit for all peaks**.



All calculated values are displayed in the Brief Report group box of the data window.



- On the Workflow button bar, click **Report**.



- In the Fitting report window click **File** and **Print** to print the report.

```

Fitting report
File Edit Search
1 1D FIT RESULTS
2 -----
3
4 Dataset : C:/p/z/data/DOSY_exp/5/pdata/1/ct1t2.txt
5
6 AREA fit : Diffusion : Variable Gradient :
7
8 I=I[0]*exp(-D*SQR(2*PI*gamma*6i*LD))*(BD-LD/3)*1e4)
9
10 25 points for Integral 1, Integral Region from 8.610 to 8.446 ppm
11
12 Converged after 49 iterations!
13
14 Results      Comp. 1
15
16 I[0]          = 9.817e-001
17 Diff Con.    = 1.585e-010 m2/s
18 Gamma       = 4.258e+003 Hz/G
19 Little Delta = 2.800m
20 Big Delta   = 99.900m
21
22 RSS         = 1.348e-005
23 SD         = 7.343e-004
24
25 Point      Gradient      Expt      Calc      Difference
26
27 1 9.630e-001  9.830e-001  9.809e-001 -2.054e-003
28 2 2.829e+000  9.742e-001  9.748e-001  6.746e-004
29 3 4.695e+000  9.623e-001  9.629e-001  5.530e-004
30 4 6.560e+000  9.446e-001  9.453e-001  6.885e-004
31 5 8.426e+000  9.222e-001  9.223e-001  6.356e-005
32 6 1.029e+001  8.942e-001  8.944e-001  1.490e-004
33 7 1.216e+001  8.619e-001  8.620e-001  1.291e-004
34 8 1.402e+001  8.254e-001  8.257e-001  3.019e-004
35 9 1.589e+001  7.862e-001  7.861e-001 -2.347e-005
36 10 1.775e+001  7.437e-001  7.439e-001  2.300e-004
37 11 1.962e+001  6.998e-001  6.996e-001 -1.402e-004
38 12 2.149e+001  6.542e-001  6.539e-001 -2.506e-004
39 13 2.335e+001  6.095e-001  6.075e-001 -1.988e-003
40 14 2.522e+001  5.610e-001  5.609e-001 -4.750e-005
41 15 2.708e+001  5.140e-001  5.148e-001  7.867e-004
42 16 2.895e+001  4.686e-001  4.695e-001  9.776e-004
43 17 3.082e+001  4.250e-001  4.257e-001  6.155e-004
44 18 3.268e+001  3.835e-001  3.835e-001  2.036e-005
45 19 3.455e+001  3.428e-001  3.434e-001  6.505e-004
1 : 1

```


- In the Fitting report window click **File** and select **Close**.

### 11.3.10 Displaying the DOSY Plot

- On the Workflow button bar, click **Back**.
- On the command line, type **re 5**.
- On the menu bar, click **Analyse**.
- On the Workflow button bar, click **Dosy**.


# Proton DOSY Experiment

- From the drop-down list select **Setup Parameters**.



General		
Method	exponential	Processing method
ExpVar	Gradient	Variable parameter
Xlist	difflist	Variable parameter values file name
Nstart	0	Start of input points
Ndata	25	Number of input points (TD)
Maxiter	100	Maximum number of iterations
EPS	1	Tolerance
Nexp	1	Number of components to fit
Noise	0	Noise level (S_DEV)
PC	4	Noise sensitivity factor
SpiSup	1	Spike suppression factor
F1mode	Peaks	F1 output data mode
Imode	Integral	Fitted intensity meaning
Scale	Logarithmic	Scaling
LWF	1	Line width factor
DISPmin	-10	Lower display limit
DISPmax	-8.02228	Upper display limit
Npars	7	Number of parameters

- Make the following change:  
Scale = **Logarithmic**
- On the toolbar, click **Start fitting**.



# 12 Multiplet Analysis

## 12.1 Introduction

---

This analysis tool can be used to define multiplets and deduce chemical shifts, coupling constants, multiplicities and connections.

## 12.2 Sample

---

100 mg 2, 3,-Dibromopropionic acid in CDCl<sub>3</sub>

## 12.3 Multiplet Assignments

---

### 12.3.1 Preparation Experiment


---


Run a 1D Proton spectrum, following the instructions in the TopSpin Guide Book *Basic NMR Experiments*, chapter *1D Proton Experiment, Experiment Setup through Processing*.

### 12.3.2 Limit Settings

---

Changing the Sweep width to a smaller value will increase the resolution.

- On the command line, type:  
**wrpa 2**  
**re 2**
- Expand the spectrum from **10.6 ppm** to **-0.5 ppm**.
- On the toolbar, click **Set sw to current region** and **SFO1 to center of region**. 




### 12.3.3 Acquisition

---


- To adjust the receiver gain, click **Gain**.
- To start the acquisition, click **Run**.

## 12.3.4 Processing


- On the Workflow button bar, click **Process**.
- On the **Proc Spectrum** button, click the **drop-down** arrow to see more options.
- In the list, select **Configure Standard Processing (proc1d)**.





- Enable the following options:
  - Exponential Multiplay (em)
  - Auto - Phasing (apk)
  - Set Spectrum Reference (sref)
  - Auto Baseline Correction (abs)
- Change LB [Hz] = 0



- In the proc1d window, click **Execute**.




- On the toolbar, click **Exact Zoom**. 
- In the exactzoom window enter the following parameters:  
From = **4.6**  
To = **3.6**
- In the exactzoom window, click **OK**.




- Right-click the spectrum window and on the shortcut menu select **Save Display Region To**.



- Enable the option **Parameters F1/2 [dp1]**.
- In the Save display region window, click **OK**.






## 12.3.5 Peak Picking

- On the menu bar, click **Analyse**.
- On the Workflow button bar, click **Pick Peaks**.


This enters the manual peak picking mode. The Dataset tabs are replaced by the Peak Picking Tool bar.



By default, the **Define new peak picking range** button is enabled.

- Click and draw a rectangle over all multiplets up to 3.7 ppm.
- On the Peak Picking toolbar, click **Modify existing peak picking range** to manually adjust the minimum and maximum intensity levels. 

- Adjust the bottom line of the box to be above the baseline (Minimum intensity) and the top line above the highest peak of all multiplets (Maximum intensity).




- On the Peak Picking toolbar, click **Return, save region** to store the peak values.




## 12.3.6 Assigning the Multiplets

- Expand the multiplet at 4.5 ppm.
- On the menu bar, click **Analyse**.
- On the **Multiplets** button, click the **drop-down** arrow to see more options.
- In the list, select **Enter multiplet analysis**.

This enters the multiplet analysis mode. The Dataset tabs bar is replaced by the Multiplet analysis button bar.



- Click **Define Multiplets Manually**. 
- Place the cursor lines to the left of the first peak of the multiplet.


### Multiplet Analysis

There will be 2 cursor lines displayed. Use the right cursor line to select a peak.

**Step 1:** Move the cursor lines slowly towards the first peak.

## Multiplet Analysis


**Step 2:** The right cursor line will stop when it gets into the center of the peak.



**Step 3:** Click left.

**Step 4:** Move the cursor lines slowly towards the second peak.

**Step 5:** The right cursor line will stop when it gets into the center of the peak.






**Step 6:** Click left.

A small marker is placed above the top of the two peaks.

**Step 7:** Move the cursor lines into the center of the two marked peaks.




**Step 8:** Right-click to open the shortcut menu.


- Define Multiplet
- Finish Current Mode
- Automatically Define Multiplet
- Automatically Define Multiplet By Region
- Define Multiplet By Region
- Define Multiplet Manually
- Define Multiplet By Free Grid
- Couple Existing Multiplets
- Define Multiplet By Coupled Grid

**Step 9:** In the list, select **Define Multiplet**.

## Multiplet Analysis




- Repeat steps 1 - 9 starting with the third peak and ending with the fourth peak.



**Step 10:** Select the Couple Existing Multiplets button.

**Step 11:** Move the cursor line into the center of the first two peaks marked 1.

**Step 12:** Click left.



**Step 13:** Move the cursor line into the center of the second two lines marked 2.





While executing the next 2 steps, the colors of the brackets over the peaks 1 and 2 change from black to red.

**Step 14:** Click left.

**Step 15:** Move the cursor into the center of the displayed multiplet.


# Multiplet Analysis



**Step 16:** Right-click to open the shortcut menu.




**Step 17:** Select **Couple Existing Multiplets**.




**Step 18:** Right-click in the spectrum window to open the shortcut menu.

**Step 19:** Select **Finish Current Mode**.




# Multiplet Analysis


- Expand the multiplet at **3.9 ppm**.
- Repeat steps 1 - 19 for this multiplet.




- Expand the multiplet at **3.7 ppm**.
- Repeat steps 1-19 for this multiplet.



- Display all 3 multiplets.




- Click **Show Multiplet Report**. 
- In the Assign Connections / Report window, click **Find Connections**.

ID	Shift [p...]	J [Hz]	M	Connection
1	4.4950	11.2575	2	J(1, 0)
		4.3790	2	J(1, 0)
2	3.9142	11.2530	2	J(2, 0)
		10.0905	2	J(2, 0)
3	3.6970	10.0977	2	J(3, 0)
		4.3790	2	J(3, 0)

Buttons: Ok, Print, Copy, Save..., Start editor, JMR, JPF, Find Connections


# Multiplet Analysis

- In the Multiplet Connection options window, click **OK**.



The connections are now assigned, and the report can be printed.

- In the Assign Connections / Report window, click **OK**.



- Click **Return, save multiplets [sret]**.






# 13 Adding a New Solvent


## 13.1 Introduction

This chapter describes the procedure how to add a new solvent to the solvent list. As an example, the solvent  $C_6D_{12}$ , Cyclohexane-d12 is used.



## 13.2 Adding Cyclohexane-d12 to the Solvent List

- At the command prompt, type **edsolv**.
- In the Edsolv window, select the **Solvents** tab.



- Right-click on the  $C_6D_6$  solvent and on the shortcut menu, select **Add new solvent**.




## Adding a New Solvent

- In the Password request window, enter the password and click **OK**.




- Add the following Solvent parameters:  
Solvent name =  $C_6D_{12}$   
Solvent description = **Cyclohexane-d12**  
Enable **Lock Solvent**  
Melting Point [K] = **279**  
Boiling Point [K] = **354**


- In the Create solvent window, click **OK**.



- In the Edsolv window, select the **Lock** tab.



- Right-click on the C<sub>6</sub>D<sub>12</sub> solvent and on the shortcut menu, select **Edit lock parameters**.



- Add the following Lock parameters:

Probe description = **Default probe**

Lock power = **-30**

Loop gain = **-10**

Loop time = **0.4**

Loop filter = **200**

Lock phase = **-1**


Lock power instep = **10**

Temperature lock power = **-20**

Shift [ppm] = **1.38**

Relative intensity = **1**

Type = **Lock**




- In the Edit lock parameters window, click **OK**.

## Adding a New Solvent


The Lock parameters will be stored for the current probe and for the selected probe only. This current probe (as defined in edhead) will then also appear in the Lock parameters list. Parameters for the other probes will remain unchanged and remain visible as probe type Generic.



- In the Edsolv window, select the **Spectrum Reference** tab.



- In the Nucleus list, select **1H**.



- Right-click on the  $C_6D_{12}$  solvent and on the shortcut menu, select **Edit spectrum reference parameters**.



- In the Edit spectrum reference parameters window, add the following Spectrum reference parameters for  $C_6D_{12}$ :  
Noise factor = **10**  
Reference shift [ppm] = **0**

Reference shift correction [ppm] = 1

Search Width [ppm] = 5

- Click **Add signal regions** and add the following regions:  
 Region 1, Lower limit [ppm] **15**, Upper limit [ppm] **1.42**  
 Region 2, Lower limit [ppm] **1.35**, Upper limit [ppm] **0.2**  
 Region 3, Lower limit [ppm] **-0.2**, Upper limit [ppm] **-3**


Lower limit [ppm]	Upper limit [ppm]	Description	Delete
15	1.42	Region 1	<input type="checkbox"/>
1.35	0.2	Region 2	<input type="checkbox"/>
-0.2	-3	Region 3	<input type="checkbox"/>

- In the Edit spectrum reference parameters window, click **OK**.

Solvent	Noise	Reference Shift [ppm]	Reference Shift Corr. [ppm]	Search Width [ppm]	Signal Regions [ppm]
Acetic	10	0	0	5	1000/180, 176/21, 19/1, -1/-1000
Acetone	10	0.9	0	5	1000/207, 205/30.8, 28.8/1, -1/-1000
<b>C6D12</b>	10	0	1	5	200/27.5, 25.5/1, -1/-5
C6D6	10	0.22	0	5	1000/129, 127/1, -1/-1000

- In the Edsol window, select Nucleus = **13C**.


## Adding a New Solvent




- Right-click on the C<sub>6</sub>D<sub>12</sub> solvent and on the shortcut menu, select **Edit spectrum reference parameters**.




- Add the following Spectrum reference parameters for C<sub>6</sub>D<sub>12</sub>  
Noise factor = **10**  
Reference shift [ppm] = **0**  
Reference shift correction [ppm] = **1**  
Search Width [ppm] = **5**



- Click **Add signal regions** and add the following regions:  
Region 1, Lower limit [ppm] **200**, Upper limit [ppm] **27.5**  
Region 2, Lower limit [ppm] **25.5**, Upper limit [ppm] **1**  
Region 3, Lower limit [ppm] **-1**, Upper limit [ppm] **-5**



- In the Edit spectrum reference parameters window, click **OK**.
- In the Edsolv window, click **Close**.



## 13.3 TopShim Solvent Parameters

Follow the instructions in the *TopShim Automatic Shimming Users Manual*, chapter *Solvents* to configure the TopShim solvent parameters.





# 14 Contact

## Manufacturer

Bruker BioSpin GmbH  
Silberstreifen 4  
D-76287 Rheinstetten  
Germany

E-Mail: [nmr-support@bruker.com](mailto:nmr-support@bruker.com)

<http://www.bruker.com>

WEEE DE43181702

## Bruker BioSpin Hotlines

Contact our Bruker BioSpin service centers.

Bruker BioSpin provides dedicated hotlines and service centers, so that our specialists can respond as quickly as possible to all your service requests, applications questions, software or technical needs.

Please select the service center or hotline you wish to contact from our list available at:

<https://www.bruker.com/service/information-communication/helpdesk.html>





 **Bruker Corporation**

[info@bruker.com](mailto:info@bruker.com)  
[www.bruker.com](http://www.bruker.com)



Order No: H147756