

TOPSPIN - data processing

- o-□ nmrafd open Folder to see file-names
 - o-□ 2103yourname open Folder to see all your measurements
- drag file in display

1H processing

- ef transformation
- ap automatic phasing

phasing-manually



put cursor in "0" or "1" depress mouse and adjust phase
save.↓

expand reference-peak (drag to expand)




indicate peak (left click) enter ppm-value



bn baseline correction is better for integrals

integrate-manually



click green to integrate, (click  blanc to expand)

drag to integrate (red lines)

if the slope is not good: \int_b interactive bias correction.

integral-calibration right-click in an integral - Calibrate current integral -enter value -ok
save.↓

Peak picking:


ppf gives peaks of Full spectrum
see if peaks of interest are shown,

pps give peaks of Displayed part of spectrum

>mi minimum intensity (standard value =0.2)

If too many peaks: mi higher (e.g. 0.5) and again ppf

If small peak of interest is not shown: mi lower (e.g. 0.01) and again ppf
(to get small peaks cy=100 , mi=0.01 , pc=0.5)



easy define ppm range

xw plot

click green dots and indicate (click) spectrum

Edit -Xmin, Xmax to enter ppm-range, (e.g. 8, -0.1), ok

1D,2D-Edit

◆ adjust amplitude

click Expand draw red box for expanding range

File -print -print

close -Plot-editor

ok close and discard all changes

Cosy processing

Put the cosy file in Topspin-display

xfb 2D-transformation (or: x (xfb , levcalc))
sym symmetrization

right-click in the top-projection – external projection setproj set projections
check or change expno

ok
right-click in the left projection – external projection
check expno
ok

referencing by use of sr

re 1 (read expno 1) (if this spectrum is calibrated)
sr ↵ ctl-c (copy: this value) enter
re 2 (read cosy expno 2)
sr ↵ ctl-v tab ctl-v (paste: this value, tab , paste this value again)

xw plot
expand draw red box to expand
green dots (click) indicate (click) spectrum

1D/2D-Edit
upper symbols for 2D

◆ adjust amplitude of 2D
lower symbols is for Top or Left projection

◆ adjust amplitude of top/left-projection


close-1D/2D-edit


File print print

Close-plot-editor , **Ok** (close and discard all changes)

13C processing

ef Fourier transformation
ap automatic phasing

manual phasing click 
put cursor in "0" or "1" and move mouse to adjust phase
save ↵

 expand reference peak
indicate peak left click enter ppm value e.g. CDCl3=77.16 ppm
for D2O sr= -4.35

ppf peaklist
see if peaklist includes smallest-peak-of-interest
if there are too many peaks >mi put in higher value , (see scale for best value)
>ppf new peaklist



easy define ppm range

xw

green dots select spectrum

1D/2D-Edit

↑ click to fit top inside scale

↓ click “ bottom “

or change amplitude

File -print -print

close Plot-editor

ok close and discard all changes

Noesy - Processing

xfb (2D-transformation)

(or: x

(xfb , levcalc))

setproj

check projections , right-click in top-projection, external projection , check or change expno , ok

, right-click in left-projection , external projection , check or change expno , ok

referencing by use of sr


re 1 (read expno 1) (if this spectrum is calibrated)

sr ↵ ctl-c (copy) enter

re 3 (read Noesy expno 3)

sr ↵ ctl-v tab ctl-v (paste , tab , paste)

Phasing

click 

select peaks , right-click , add

click **R** to phase rows

adjust phasing with **0** and **1**

save ↵

click **C** to phase columns

adjust phasing with **0** and **1**

save ↵

referencing and plot same as cosy.

HSQC - Processing

xfb (2D-transformation) (or: x (xfb , levcalc))
(setproj)

check projections , right-click in **top**-projection, external projection , check or change **expno 1H**, ok
, right-click in **left**-projection , external projection , check or change **expno 13C**, ok

referencing by use of **sr** values

re 1 (read **1H**-spectrum expno 1) (spectrum must be calibrated)

sr ↵ ctl-c (copy) , enter

re 3 (read hsqc expno 3)

sr ↵ ctl-v (paste) , ↵ sr (**F2** , F1)


re 2 (read **13C**-spectrum expno 2) (spectrum must be calibrated)

sr ↵ ctl-c (copy) , enter

re 3 (read hsqc expno 3)

sr ↵ tab , (move to F1 position) , ctl-v , ↵ sr(**F2** , **F1**)

Phasing

click 

select peaks , right-click , add

click **R** to phase rows

adjust phasing with **0** and **1**

save

click **C** to phase columns

adjust phasing with **0** and **1**

save

↵

xw xwinplot

expand draw red box to expand

green dots (click) indicate (click) spectrum

1D/2D-Edit

upper symbols for 2D

◆ adjust amplitude of 2D

lower symbols is for Top or Left projection

◆ adjust amplitude of top/left-projection

close-1D/2D-edit

File print print

Close-plot-editor

Ok (close and discard all changes)